

The manual provided is designed to cover the standard configurations for the chamber you ordered. Any special features or optional items ordered are presented below.

Your chamber is equipped with the following:

- Temperature ranges of -10°C to 44°C \pm 0.5°C (lights off) and 0°C to 44°C \pm 0.5°C (lights on)
- Humidity range up to $85\% \pm 10\%$ over 15° C to 30° C (lights off).
- IntellusUltra Controller Tablet / Touch Screen App
- Light intensity of 405 μmoles / m² / sec as measured
 6" from the light fixture at a temperature of 24°C and a light lifetime of 100 hours
- SciWhiteTM LED tiles *
- Open loop dimmable lighting *
- Extended temperature range *
- Pan-type humidifier with Vaisala HMP60 relative humidity sensor
- Unit cooler drain line routed out rear wall of chamber *
- Air-cooled condensing unit

Your chamber does not include the following options as described in the manual:

- No remote alarm contacts are provided
- No dehumidification system is provided

Note: All supplemental information supersedes any similar information provided in the standard manual.