


REFRIGERATION SYSTEM CHARACTERISTICS		
VOLTAGE:	220 VAC	PHASE: 1Ø
		FREQUENCY: 50 Hz
MODEL	REFRIGERANT	CHARGE
30 SERIES	R-134A	10 oz.
30 SERIES	R-513A	9 oz.


PERCIVAL SCIENTIFIC, INC.
 TITLE: REFRIGERATION DIAGRAM - GENERAL
 MODEL: 30 SERIES DRAWING: 032156-1-1 REV: -
 DATE: 1/23/2022