

## **TEST CERTIFICATION**

This is to certify that the <u>I-66NL</u> chamber serial number <u>32801.14.23</u> provided by Percival Scientific, Inc. for Percival Scientific Hong Kong Co., Ltd. has successfully tested to meet the specifications listed below:

✓ The chamber has met the temperature range test of  $2^{\circ}$ C to  $44^{\circ}$ C ± 0.5°C \*

✓ The chamber has met the humidity range test of down to  $30\% \pm 5\%$  @  $4^{\circ}C^{*}$ 

\* - These temperatures and humidity levels are achievable when ambient conditions are 24°C  $(75^{\circ}F)$  or below. Higher ambient temperatures may affect the operating ranges.

Testing was performed with a calibrated circular chart recorder.

Date \_\_\_\_\_11/06/23 \_\_\_\_