



## TEST CERTIFICATION

This is to certify that the I-30L chamber serial number 32838.08.23 provided by Percival Scientific, Inc. for Oriental Science & Technology Co. (OSTC) has successfully tested to meet the specifications listed below:

✓ The chamber has met the temperature range test of  $2^{\circ}\text{C}$  to  $44^{\circ}\text{C} \pm 0.5^{\circ}\text{C}$  (lights off) \*  
 $7^{\circ}\text{C}$  to  $44^{\circ}\text{C} \pm 0.5^{\circ}\text{C}$  (lights on) \*

✓ The chamber has met the humidity range test of  $50\%$  to  $85\% \pm 10\%$  (lights off) \*

\* - These temperatures and humidity levels are achievable when ambient conditions are  $24^{\circ}\text{C}$  ( $75^{\circ}\text{F}$ ) or below. Higher ambient temperatures may affect the operating ranges.

Testing was performed with a calibrated circular chart recorder.

Date 10/25/2023