

Partial Parts List For Job #: 21077-0000

Date Printed: December 5, 2023

For Sales Order: 32994 Item: ARC36L4CUSWC8

| Part # | Description | QTY | UOM |
|---------------|--|------|-----|
| 2015-1C-30 | 36 SERIES INSULATED BOX (PURCHASED) REV. B IS THE NEW BOX DESIGN OFF THE SALVAGNINI WITH NEW FLOOR DRAIN & BREAKER REDESIGN. CHANGE ON BOXES SHIPPED FROM LEER AFTER 9-1-16. | 1.00 | EA |
| 2017-5-230-1 | SHORT (NORMAL) SENSOR BRACKET Bend Location: B-- | 1.00 | EA |
| 2018-8-600-1 | DOUBLE FLARE COPPER TUBE 1/2" OD COPPER, 9" LENGTH | 1.00 | EA |
| 2020-11-120-1 | 36 LAMP BANK BASE Bend Location: A-- | 4.00 | EA |
| 2020-11-120-2 | 36 LAMP BANK COVER Bend Location: A-- | 4.00 | EA |
| 2020-11-300-1 | UPPER LAMP BANK FIXTURE SUPPORT Bend Location: B-- | 2.00 | EA |
| 2020-11-70-2 | UNIT COOLER AIR BLOCK Bend Location: B-- | 1.00 | EA |
| 2020-11-70-3L | LH UNIT COOLER AIR BLOCK HANGER Bend Location: B-- | 1.00 | EA |
| 2020-11-70-3R | RH UNIT COOLER AIR BLOCK HANGER Bend Location: B-- | 1.00 | EA |
| 2020-11-70-4 | AIR BLOCK BACKING STRIP 22 GA #304-4 S.S. | 2.00 | EA |
| 2020-11-70-5 | UNIT COOLER AIR BLOCK RUBBER 1/16" X 4 3/4" X 6" NEOPRENE RUBBER | 2.00 | EA |
| 2020-11-90-1 | FIXTURE SUPPORT Bend Location: B-- | 6.00 | EA |
| 2020-11-90-2 | AIR PLENUM SUPPORT Bend Location: B-- | 6.00 | EA |
| 2020-11-90-2 | AIR PLENUM SUPPORT Bend Location: B-- | 2.00 | EA |
| 2020-11-95-1 | BOTTOM FIXTURE SUPPORT 22 GA #304-4 S.S. | 2.00 | EA |
| 2020-2-100-2 | 36 CONTROL BOX ENCLOSURE Bend Location: A-- | 1.00 | EA |
| 2020-2-100-3 | 36 CONTROL BOX TOP COVER Bend Location: B-- | 1.00 | EA |
| 2020-2-100-4 | 36 CONTROL BOX DOOR STOP Bend Location: B-- | 1.00 | EA |
| 2020-2-100-5 | 36 SOLID-STATE RELAY MOUNTING PLATE Bend Location: B-- | 1.00 | EA |
| 2020-2-100B | CONTROL BOX LOCK CATCH PLATE WELD | 1.00 | EA |
| 2020-2-100C | 36 CONTROL BOX DOOR WELD ASS'Y | 1.00 | EA |
| 2020-2-80-1L | LH CONTROL TOP CORNER Bend Location: B-- | 1.00 | EA |
| 2020-2-80-1R | RH CONTROL TOP CORNER 18 GA PTLOK | 1.00 | EA |
| 2020-2-80-2 | 36 CONTROL TOP REAR PANEL 22 GA PTLOK | 1.00 | EA |
| 2020-2-80-3 | 36 MECH SECTION COVER 22 GA PTLOK | 1.00 | EA |
| 2020-21-110-1 | DRIVER RACK Bend Location: B-- | 1.00 | EA |
| 2021-8-170-1 | HEATER BRACKET | 2.00 | EA |
| 57458 | DRIER CLAMP Bend Location: B-- | 1.00 | EA |
| 70-2-300-16 | SUCTION ACCUMULATOR MTG. BRACKET Bend Location: B-- | 2.00 | EA |
| 93-5-140-1 | AIR PLENUM Bend Location: A-- | 3.00 | EA |
| 93-5-140-1 | AIR PLENUM Bend Location: A-- | 1.00 | EA |
| 93-5-140-2 | AIR DIFFUSER PLATE Bend Location: A-- | 3.00 | EA |
| 93-5-140-2 | AIR DIFFUSER PLATE Bend Location: A-- | 1.00 | EA |
| 93-5-140-5 | AIR DEFLECTOR 23% OPEN ALUM. Bend Location: A-- | 3.00 | EA |
| 93-5-140-5 | AIR DEFLECTOR 23% OPEN ALUM. Bend Location: A-- | 1.00 | EA |
| 95-3-120-3 | SIDE PANEL CHANNEL Bend Location: B-- | 2.00 | EA |
| 95-3-120-3 | SIDE PANEL CHANNEL Bend Location: B-- | 2.00 | EA |
| 95-3-300-2 | CONDENSING UNIT PLATE Bend Location: B-- | 1.00 | EA |
| 95-3-300-5 | TOP LH PANEL Bend Location: B-- | 1.00 | EA |
| 95-3-300-6 | TOP RH SCREEN Bend Location: B-- | 1.00 | EA |
| 95-4-140-3 | AIR DIFFUSER SCREEN Bend Location: A-- | 3.00 | EA |
| 95-4-140-3 | AIR DIFFUSER SCREEN Bend Location: A-- | 1.00 | EA |
| 95-4-140-4 | LOUVER | 3.00 | EA |

Partial Parts List For Job #: 21077-0000

Date Printed: December 5, 2023

For Sales Order: Item: ARC36L4CUSWC8

| Part # | Description | QTY | UOM |
|----------------------|--|-------|-----|
| 95-4-140-5 | <i>Bend Location: B--</i> BOTTOM LOUVER W/ NOTCH 22 GA PTLOK | 1.00 | EA |
| 95-4-300-3 | REAR AIR BLOCKS 22 GA PTLOK 11 15/16 X 26 | 4.00 | EA |
| 95-4-300-4 | SHIPPING CLIP-AIR PLENUM 22 GA. #304-5 S.S. 3/4" X 2 15/16" | 2.00 | EA |
| 95-4-300-97 | RH RECEPTACLE PANEL <i>Bend location: A--</i> | 1.00 | EA |
| 95-4-50-1 | DOOR FACE PANEL <i>Bend Location: A--</i> | 1.00 | EA |
| 95-4-50-2 | DOOR CORE PANEL <i>Bend Location: A--</i> | 1.00 | EA |
| 95-4-50-5 | CORNER-DOOR CORE <i>Bend Location: FLAT</i> | 4.00 | EA |
| 95481 | SPORLAN SOLENOID VALVE MOUNT 16 GA GALV 1 5/8" X 3 3/8" | 2.00 | EA |
| 98-2-300-7 | LH RECEPTACLE PANEL 22 GA PTLOK | 1.00 | EA |
| CMS-027 (MACH. | 1/4-20 UNC X 1/2", ROUND HD., ZINC ZINC-CHROMATE PLATED | 6.00 | EA |
| CMS-044 (MACH. | 6-32 X 1/2" PHILLIPS ROUND HD. S.S. | 4.00 | EA |
| CMS-044 (MACH. | 6-32 X 1/2" PHILLIPS ROUND HD. S.S. | 4.00 | EA |
| CMS-044 (MACH. | 6-32 X 1/2" PHILLIPS ROUND HD. S.S. | 3.00 | EA |
| CMS-069 (MACH. | 10-24 UNC X 3/8, PH. RD. HD., S.S. 10-24 UNC X 3/8" PHILLIPS, ROUND HEAD, STAINLESS-STEEL MACHINE SCREW | 2.00 | EA |
| CMS-074 (MACH. | 6-32 X 1" PHILLIPS RD. HD., S.S. | 4.00 | EA |
| CMS-111 (MACH. | 6-32 X 1/4" PHIL. PAN HEAD, ZINC LED LIGHT MOUNTING SCREWS: Thrd-Cutting Screw for Metal&Plastic (Type F) Pan Head Phil, Zinc-Plated Steel, 6-32 Thread, 1/4" L | 64.00 | EA |
| CSM-004 (S.M. SCREW) | 6 X 1/2" NICKEL PL, PHIL. TRUSSHEAD #6128-148N SERVICE SUPPLY CO. | 4.00 | EA |
| CSM-004 (S.M. SCREW) | 6 X 1/2" NICKEL PL, PHIL. TRUSSHEAD #6128-148N SERVICE SUPPLY CO. | 4.00 | EA |
| CSM-020 (S.M. SCREW) | 10-16 X 3/4" HEX WASHER HEAD TEKS S 10-16 X 3/4" HEX WASHER HEAD SHEET METAL SCREW, DRILL POINT, GRADE 2, ZINC-PLATED STEEL | 24.00 | EA |
| CSM-039 (S.M. SCREW) | 6 X 1/2" SELF DRILLING #410 S.S. HEX HEAD, CORROSION RESISTANT COATING | 15.00 | EA |
| CSM-039 (S.M. SCREW) | 6 X 1/2" SELF DRILLING #410 S.S. HEX HEAD, CORROSION RESISTANT COATING | 36.00 | EA |
| CSM-039 (S.M. SCREW) | 6 X 1/2" SELF DRILLING #410 S.S. HEX HEAD, CORROSION RESISTANT COATING | 40.00 | EA |
| CWS-024 (DRYWALL ... | 8 X 2 PARTICAL BOARD SCREW PARTIAL THREAD 8 X 2 PHILIPS HEAD FLAT HEAD WITH NIBS, BLACK PHOS STEEL | 41.00 | EA |
| EBR-145 (BREAKER) | 10A, 1-P, ALTECH 1DU10L, UL489 ALTECH CORP. CATALOG #1DU10L, DIN MOUNT (AC) CIRCUIT BREAKER, 10A, 1-POLE, 10kA @ 120/240/277VAC, 50/60Hz, L-SERIES; SCREW TERMINALS, GREEN: OFF/TRIP, RED: ON, LINE/LOAD REVERSIBLE, HIGH MAGNETIC TRIP POINT, UL 489 & CE | 1.00 | EA |
| EBR-146 (BREAKER) | 15A, 1-P, ALTECH 1DU15L, UL489 ALTECH CORP. CATALOG #1DU15L, DIN MOUNT (AC) CIRCUIT BREAKER, 15A, 1-POLE, 10kA @ 120/240/277VAC, 50/60Hz, L-SERIES; SCREW TERMINALS, GREEN: OFF/TRIP, RED: ON, LINE/LOAD REVERSIBLE, HIGH MAGNETIC TRIP POINT, UL 489 & CE | 1.00 | EA |
| EBR-147 (BREAKER) | 2A, 1-P, ALTECH 1DU2L, UL489 ALTECH CORP. CATALOG #1DU2L, DIN MOUNT (AC) CIRCUIT BREAKER, 2A, 1-POLE, 10kA @ 120/240/277VAC, 50/60Hz, L-SERIES; SCREW TERMINALS, GREEN: OFF/TRIP, RED: ON, LINE/LOAD REVERSIBLE, HIGH MAGNETIC TRIP POINT, UL 489 & CE | 4.00 | EA |
| EBU-016 (BUSHING) | SR7N-2 HEYCO INSULATOR BUSHING HEYCO #1244 STRAIN RELIEF BUSHING; BLACK, UL/CSA STRAIGHT THRU-TYPE; FOR USE W/ Ø 0.430" MAX. ROUND CABLES; FOR USE W/ (EPL-030) 7-WIRE CORDSET. REQUIRES A 0.66" WIDE X 0.75" LONG OVAL MOUNTING HOLE. | 4.00 | EA |
| ECO-010 (CONNECTOR) | 3/8" STRAIGHT ROMEX 2711 RACO 2800 I.E. 6623, NON-METALLIC CABLE, REQUIRES ONE (1) 7/8" DIA. MOUNTING HOLE; UL- & CSA-APPROVED; KRIZ. STCNC301 | 1.00 | EA |
| ECO-061 (CONNECTOR) | WAGO 231-304/026-000 FOUR CIRCUIT WAGO #231-304/026-000; 4-POLE ORANGE FEMALE CONNECTOR W/ CAGE CLAMP, FOR USE W/ AWG #28-12 WIRE, 300 V/15 A MAX., UL & CSA APPROVED. | 2.00 | EA |
| ECO-119 (CONNECTOR) | WAGO 2 CONDUCTORS LEVER CONNECTOR WAGO 221-412, 24-12 GA. WIRE, 600V, 20A, 2 COMMON TERMINALS, UL, CE, EN, 100/BOX | 18.00 | EA |
| ECO-119 (CONNECTOR) | WAGO 2 CONDUCTORS LEVER CONNECTOR WAGO 221-412, 24-12 GA. WIRE, 600V, 20A, 2 COMMON TERMINALS, UL, CE, EN, 100/BOX | 12.00 | EA |
| ECO-120 (CONNECTOR) | WAGO 3 CONDUCTORS LEVER CONNECTOR WAGO 221-413, 24-12 GA. WIRE, 600V, 20A, 3 COMMON TERMINALS, UL, CE, EN, 50/BOX | 1.00 | EA |
| ECO-121 (CONNECTOR) | WAGO 5 CONDUCTORS LEVER CONNECTOR WAGO 221-415, 24-12 GA. WIRE, 600V, 20A, 5 COMMON TERMINALS, UL, CE, EN, 25/BOX | 5.00 | EA |
| EHE-046 (HEATER) | 250W @ 230V U-TUBE HEATER, RED OPEN AIR HEATER, 250W @ 230V, MARKED WITH RED PAINT, WATLOW PART RBN36B1H-2, 40" OAL X 0.315" DIA., FORMED IN FIGURE 1, A = 17.5625", R = 1", W/ EPOXY SEAL AND 36" SILICONE WIRE LEADS W/ SHRINK BOOT EPOXY SEAL, QUOTE # 710091 See drawing number 2006-1-300-7 | 1.00 | EA |
| EIN-002 (INSULINER) | 422 T&B 1/2", BLUE OR BLACK THOMAS & BETTS #422; 1/2" TRADE SIZE, BLUE OR BLACK, INSULINER; UL-APPROVED (UL FILE #E23018); FOR USE W/ (ECO-010). | 2.00 | EA |
| ELD-049 | DALI TO ANG CONVERTER (0-10VDC OUT) LT-84A DIMMING SIGNAL CONVERTER; DALI, PUSH DIM SIGNAL INPUT, 0-10V ANALOG x 4CH | 1.00 | EA |
| ELD-066 | 12 X 24 WHITE LED TILE W/COATING 12" x 24" LED tile with conformal coating. With Samsung LM301B Mid Power LED, 4000K with CR190 SH flux bin. 12 STRINGS OF LEADS WITH 13 LEADS PER STRING (156). | 8.00 | EA |
| ELD-067 | LED DRIVER, MEANWELL HLG-60H-C700B | 4.00 | EA |

Partial Parts List For Job #: 21077-0000

Date Printed: December 5, 2023

For Sales Order: Item: ARC36L4CUSWC8

| Part # | Description | QTY | UOM |
|------------------------|--|--------|-----|
| | MEANWELL HLG-C700B LED DRIVER, CONSTANT CURRENT, 700mA OUTPUT CURRENT, 50-100V OUTPUT, 90-305VAC INPUT VOLTAGE, WITH DIMMING, -40 TO 80°C OPERATING TEMP RANGE, RoHS, 5.80"L x 2.42"W x 1.42"H, IP65 RATING, 3 IN 1 DIMMING | | |
| EPC-032 (USB CABLE) | PANEL MOUNT USB CABLE, A-A F/M, 3FT STARTECH USBPNLAFAM3 | 1.00 | EA |
| EPL-006 (PLUG) | MALE PLUG W/LEAD 18/3 W/54"LD WHITE EMF CORP. M-3-7/PM0030; 3-WIRE CORDSET, W/ 54" LENGTH AWG #18 STRANDED Cu WIRES & 300 V/105°C MAX., 1/2" PRESTRIPPED WHITE PVC INSULATION; UL- & CSA-RECOGNIZED. | 4.00 | EA |
| EPL-041 (PLUG) | 9 TERM., AMP MATE-N-LOK 1-480706-0 TYCO UNIVERSAL MATE-N-LOK 1-480706-0; 9-TERMINAL PLUG, 19 A @ 600 V MAX., 105°C MAX.; USE W/ (ERE-046); UL- & CSA-APPROVED. | 3.00 | EA |
| EPL-059 (PLUG) | TYCO 172167-1 MINI M-N-L , 4 SOCKET 172167-1 MINI MATE-N-LOK Connector, Soft Shell; 4; 600 VAC; Plug; Double-Row; Nylon; Plug; Natural | 1.00 | EA |
| EPS-031 (POWER SU...) | 120/240V IN, 12VDC OUT, 75W MEAN WELL LRS-75-12 | 1.00 | EA |
| ERE-014 (RECEPTACLE) | FEMALE PLUG W/LEAD, 18/3W/12"LD EMF CORP. PM0020; 3-WIRE FEMALE RECEPTACLE, W/ THREE (3) 12" LENGTH X AWG #18 STRANDED Cu WIRES & 300 V/105°C MAX., 1/2" PRESTRIPPED WHITE PVC INSULATION; REQUIRES A 1" X 1" SQUARE MOUNTING HOLE; FITS A 0.060" MAX. PANEL THICKNESS; UL- & CSA-RECOGNIZED. | 4.00 | EA |
| ERL-030 (RELAY) | DPDT,30A.,240VAC | 1.00 | EA |
| ERL-060 (RELAY) | S.S., N.O., 3-32 VDC INPUT, 5A, DIN CARLO GAVAZZI RP1A23D5M1; N.O., 3-32 VDC INPUT, 24-240 VAC OUTPUT, 5 A MAX. LOAD CURRENT, ZERO-SWITCH TRIAC; DIN MOUNT, PUSH-TO-RELEASE WIRING CONNECTIONS; UL- & CSA-APPROVED. | 3.00 | EA |
| ESW-130 (SWITCH) | DPST ROCKER SWITCH, 20A, 125V PRCDA1-20F-BB08W DPST ROCKER SWITCH, ON-OFF, 20A @ 125V, 10A @ 250V, PANEL MOUNT, SNAP-IN, BLACK, I-O MARKINGS, RoHS COMPLIANT, | 1.00 | EA |
| ETC-052 (TEMP. CON...) | 3 WIRE RTD, 100 OHMS A/100-3W-2.5"-3/16"-10'CL2P-HS, 120° LEADS, 3 WIRE RTD SENSOR | 2.00 | EA |
| ETC-108 (TEMP. CON...) | INTELLUS CONTROL OVERLAY - CHINA INTELLUS CONTROL OVERLAY WITH GREEN BACKGROUND FOR CHINA | 1.00 | EA |
| ETC-109 (TEMP. CON...) | INTELLUSULTRA GEN2, LOW END VER2 WATLOW 2042-4104 | 1.00 | EA |
| ETE-017 (TERMINAL) | WAGO 2004-1401, GRAY, 4 TERMINAL 5 mm TERMINAL BLOCK WIDTH FOR MOUNTING ON DIN RAIL | 8.00 | EA |
| ETE-018 (TERMINAL) | WAGO 2004-1404, BLUE, 4 TERMINAL 5 mm TERMINAL BLOCK WIDTH FOR MOUNTING ON DIN RAIL | 4.00 | EA |
| ETE-019 (TERMINAL) | WAGO 2004-1407, GROUND, 4 TERMINAL 5 mm TERMINAL BLOCK WIDTH FOR MOUNTING ON DIN RAIL | 2.00 | EA |
| ETE-020 (TERMINAL) | WAGO 2004-1491, END PLATE 5 mm TERMINAL BLOCK WIDTH FOR MOUNTING ON DIN RAIL | 2.00 | EA |
| ETE-021 (TERMINAL) | WAGO 2004-402; GRAY TERM. JUMPER WAGO #2004-402; LIGHT GRAY, 2-PIN PUSH-IN JUMPER BAR; FOR USE W/ 5 mm WIDE TERMINAL BLOCKS; 32 A MAX. | 4.00 | EA |
| ETE-022 (TERMINAL) | WAGO 249-116, END TERMINAL ANCHOR FOR 5 mm TERMINAL BLOCK WIDTH | 8.00 | EA |
| EWR-002 (TYRAP) | 6" CABLE TIES 100/BAG | 30.00 | EA |
| EWR-002 (TYRAP) | CABLE TIES, RP14897 | 4.00 | EA |
| HBS-001 (SPACER) | FEMALE THREADED CIRCUIT BOARD POST MICRO PLASTICS INC #27FTP00250 | 6.00 | EA |
| HHA-008 (HANDLE) | BLACK SOUTHCO P1-30-101-11 ABS LOW GLOSS | 1.00 | EA |
| HHI-029 (HINGE) | SOUTHCO EH-5A-4G4-50 SERIES 40 BLACK POWDER COATED MOUNTING METHOD: M5 FOUR THREADSTUDS AT .984" BETWEEN CENTER. | 2.00 | EA |
| HLO-019 (LOCK) | COMPRESSION LOCK (1/4 TURN) 040-1-1-1-42-59, ELECTRAUSA chrome with slotted lock key | 1.00 | EA |
| HNU-038 (NUT) | HEXNUT, 6/32 ONLY ZINC CHROMATE | 15.00 | EA |
| HNU-052 (NUT) | M5, .8 COARSE THREAD NUT | 8.00 | EA |
| HNU-068 (NUT) | LOCKNUT, 6-32, 18-8 S.S. McMASTER-CARR #91831A007 | 4.00 | EA |
| HNU-076 (RIVET NUT) | 6-32 RUBBER-INSULATED RIVET NUT McMaster-Carr 93495A120, Rubber-Insulated Brass Rivet Nut, 6-32 Thread, 1/64" - 5/32" Material Thick Range | 4.00 | EA |
| HST-002 (PORT PLUG) | 1 1/4" ACCESS PORT PLUG SHAW AGRO DEVELOPMENT #51013 | 1.00 | EA |
| HWA-008 (WASHER) | 1/4" FLAT, ZINC PLATED | 6.00 | EA |
| HWA-011 (WASHER) | 10 EXTERNAL LOCK ZINC-PLATED STEEL | 2.00 | EA |
| HWA-012 (WASHER) | 1/4" LOCK WASHER ZINC PLATED FASTENAL 1133618 | 6.00 | EA |
| LST-007 (STAPLES) | 15/16 CROWN X 3/4 LEG 16GA. STAPL ROOFING STAPLE | 144.00 | EA |
| LST-008 (STAPLES) | 15/16 CROWN X 1 LEG 16GA. STAPLE INTERCHANGE WC GALVANIZED STAPLE | 267.00 | EA |
| RAC-006 (SUCTION A...) | 3816 REF. RESEARCH, 1/2 ODF REFRIGERATION RESEARCH #3816; Ø 4" X 6-1/2" L, 1/2" ODF SOLDER INLET & OUTLET CONNECTIONS; 2.3# of R-134a HOLDING CAPACITY. | 1.00 | EA |
| RCO-052 (UNIT COOLER) | BOHN TA17SG, 208-230 VAC, 60 HZ BOHN T17SG; 208-230 VAC/1-PHASE/60 Hz; W/ 1/2 SAE FLARE NUT INLET & 3/8 ODF SOLDER SUCTION-LINE CONNECTIONS, 1/2 O.D. CONDENSATE DRIP PAN DRAIN CONNECTION; W/ 1,700 BTU/(hr.-10°F TD) COOLING CAPACITY; R-134a, W/ TWO (2) 208-230 VAC/1-PHASE/60 Hz FANS, 0.8 FLA TOTAL @ 208-230 VAC; 210 cfm TOTAL; UL & cUL LISTED; NSF 7-1997 CLASSIFIED. | 1.00 | EA |

Partial Parts List For Job #: 21077-0000

Date Printed: December 5, 2023

For Sales Order: Item: ARC36L4CUSWC8

| Part # | Description | QTY | UOM |
|--------------------------|---|-------|-----|
| RCU-313 | M7ZH-0045-IAZ-328 1/2HP 2022: 1/2 HP, R-134A, 220/240, 1 phase, 50 Hz. R.L.A= 3.9, MCA= 5.4 AMPS, 2.5 lbs Receiver @ 90% full (R-134a) | 1.00 | EA |
| RDR-008 (DRIER) | SPORLAN C-052-S DRIER 1/4" LINE | 1.00 | EA |
| RES-001 (ELBOW) | ELBOW, 3/8NOM C FTG WL1/2"OD WRO W-02321, STREET FITTING | 1.00 | EA |
| RFN-001 (FLARE NUT) | FREY #14FS4, 1/4" SHORT TYPE FORGED | 2.00 | EA |
| RFN-004 (FLARE NUT) | NST4-8 1/2"OD TUBE, STA-TITE | 1.00 | EA |
| RLI-001 (LIQ. INDICATOR) | SPORLAN SA-12S 1/4 ODF SIGHT GLASS 25 per box | 1.00 | EA |
| RPR-001 (PRES. REG...) | CROT-6 5/8"ODF W/STRAINER 0/60 CRANKCASE PRESSURE REGULATOR | 1.00 | EA |
| RTE-030 (TEE) | W-4002, 1/2" X 1/2" X 3/8" | 1.00 | EA |
| RVH-002 (HAND VALVE) | A-14850C MUELLER 3/8" ODF | 1.00 | EA |
| RVQ-001 (BODY) | Q BODY 2 X 4 SAE 168011, 1/4 SAE INLET, 1/2 SAE OUTLET | 1.00 | EA |
| RVQ-102 (CARTRIDGE) | QC-2 GREEN Q CARTRIDGE/ID TAG 168074 | 1.00 | EA |
| RVQ-200 (ELEMENT) | KT-43-JC 30' ELEMENT 180314 | 1.00 | EA |
| RVS-027 (VALVE) | FEMALE DIN CONNECTOR ASSY GOES W/RVS-028 ON RVS-030 & RVS-031 | 2.00 | EA |
| RVS-028 (VALVE) | OMKC-1E 220-240/50 COIL GOES W/RVS-027 ON RVS-030 & RVS-031 | 2.00 | EA |
| RVS-030 (VALVE) | SOLENOID VALVE #XWG-2, LESS COIL SPORLAN #XWG; W/ EXTENDED 3/8" ODF SOLDER CONNECTIONS | 1.00 | EA |
| RVS-031 (VALVE) | SOLENOID VALVE #XUP-7, LESS COIL SPORLAN #XUP; W/ EXTENDED 3/8" ODF SOLDER CONNECTIONS | 1.00 | EA |
| SWS-041 (SHELF) | SHELF, 27" X 28 3/4" W/WHITE EPOXY FINISH, DWG. 95-3-300-99EP | 4.00 | EA |
| XBU-001 (BUMPER) | 0.78" O.D. X 0.20" THICK, CLEAR BUMPER SPECIALTIES #BS-41, 0.78" O.D. X 0.20" THICK, CLEAR, SELF-ADHESIVE BACK | 8.00 | EA |
| XGA-040 (GASKET) | 31 1/2" X 59 3/4"OD, 4 SIDED BLACK MAGNETIC DOOR GASKET | 1.00 | EA |
| XGR-002 (GROMMET) | 3/4" I.D. X 7/8" HOLE | 2.00 | EA |
| XGR-002 (GROMMET) | 3/4" I.D. X 7/8" HOLE | 2.00 | EA |
| XGR-002 (GROMMET) | 3/4" I.D. X 7/8" HOLE LAVELLE INDUSTRIES #1022 | 2.00 | EA |
| XGR-003 (GROMMET) | 1-1/4 ID X 1-1/2 HOLE LAVELLE #1047 LAVELLE #1047, 1-1/4" I.D. X FOR 1-1/2" DIA. HOLE X 3/8" THICKNESS X 1/16" GROOVE WIDTH, BLACK | 2.00 | EA |
| XGR-003 (GROMMET) | 1-1/4 ID X 1-1/2 HOLE LAVELLE #1047 LAVELLE #1047, 1-7/8" O.D. X 1-1/4" I.D. X Ø 1-1/2" PANEL HOLE X 3/8" O.A. THICKNESS X 1/16" GROOVE WIDTH, BLACK RUBBER. | 1.00 | EA |
| XGR-003 (GROMMET) | 1-1/4 ID X 1-1/2 HOLE LAVELLE #1047 LAVELLE #1047, 1-7/8" O.D. X 1-1/4" I.D. X Ø 1-1/2" PANEL HOLE X 3/8" O.A. THICKNESS X 1/16" GROOVE WIDTH, BLACK RUBBER. | 1.00 | EA |
| XGR-005 (GROMMET) | 3/8" I.D. X 1/2" HOLE, BLACK VINYL BLACK VINYL, G1089 PLASTIC | 2.00 | EA |
| XGR-007 (GROMMET) | 5/8" ID X 7/8" HOLE #875 PLASTIC BLACK VINYL | 6.00 | EA |
| XGR-022 (GROMMET) | WHT LCKING GROMMET FOR 3/8 DIA HOLE McMaster-Carr 5269T232, Locking Grommet for 3/8" Hole Diameter, 3/16" ID, White | 8.00 | EA |
| XLA-013 (LABEL) | WAGO 44 X 99 MM, SILVER LABEL WAGO 210-804, 44 MM X 99 MM, 500 MARKERS PER ROLL, WEIGHS 400 GRAMS PER ROLL, SILVER-COLORED | 2.00 | EA |
| YCA-005 (CAPLUG) | BPF-1-5, 1" HOLE | 6.00 | EA |
| YDR-002 | ACCESS PORT SLEEVE, 1 1/4" ID COMPONENT HARDWARE R90-1002 Flexible Insulating Conduit Sleeve | 1.00 | EA |
| YHI-002 (HINGE) | IN LINE KNUCKLE STYLE HINGE SOUTHCO 96-10-400-11 | 3.00 | EA |
| YHS-001 (SPACER) | HINGE SPACER 95-3-50-5A | 3.00 | EA |
| YWA-001 (WASHER) | 30% GLASS FILLED NYLON/ BLACK .040 PLASTIC SHOULDER WASHER 24SW0006 0.040" Nylon 6/6 Shoulder Washer with Natural Finish, Off White | 64.00 | EA |