

The manual provided is designed to cover the standard configurations for the chamber you ordered. Any special features or optional items ordered are presented below.

Your chamber is equipped with the following:

- Temperature ranges of -15°C to $44^{\circ}\text{C} \pm 0.5^{\circ}\text{C}$ (lights off) and -10°C to $44^{\circ}\text{C} \pm 0.5^{\circ}\text{C}$ (lights on)
- IntellusUltra Controller Tablet / Touch Screen App
- Light intensity of $515 \mu\text{moles} / \text{m}^2 / \text{sec}$ as measured 6" from the light fixture at a temperature of 24°C and a light lifetime of 100 hours
- SciWhite® LED tiles *
- Open loop dimmable lighting *
- Unit cooler drain line routed out rear wall of chamber *
- Air-cooled condensing unit

Your chamber does not include the following options as described in the manual:

- No remote alarm contacts are provided
- No humidification or dehumidification system is provided.

Note: All supplemental information supersedes any similar information provided in the standard manual.