

Partial Parts List For Job #: 21244-0000

Date Printed: March 13, 2024

For Sales Order: 33283 Item: LT41VLSWC8T

Part #	Description	QTY	UOM
2006-1-20-2A	OUTER SHELL BOTTOM 18 GA PTLOK, 38 13/16" X 35 5/8"	1.00	EA
2006-1-20-3	OUTER SHELL REAR WALL 22 GA PTLOK, 59 9/16" X 38 7/8"	1.00	EA
2006-1-20-6A	OUTER SHELL TOP 18 GA PTLOK, 38 13/16" X 35 5/8"	1.00	EA
2006-1-20-8	TOP / BOTTOM LINER SUPPORT Bend Location: B--	2.00	EA
2006-1-20-9	VERTICAL LINER SUPPORT Bend Location: A--	2.00	EA
2006-1-30-1	INNER LINER BOTTOM 22 GA, #304-4 S.S., 38 5/16" X 29 5/32"	1.00	EA
2006-1-30-2	INNER LINER TOP Bend Location: B--	1.00	EA
2006-1-30-5	INNER LINER BACK Bend Location: A--	1.00	EA
2006-1-30-6	41 TOP & BOTTOM INSULATION POLYSTYRENE 1 9/16 x 14 7/8 x 37 13/16	4.00	EA
2006-1-50-1	DOOR FACE PANEL Bend Location: A--	1.00	EA
2006-1-50-3	CU-41 DOOR INSULATION POLYSTYRENE, 1 7/8 X 39 X 59 1/2"	1.00	EA
2006-1Q-20	BASE ASSY	1.00	EA
2006-4-20-1L	L.H. SIDE WALL Bend Location: A--	1.00	EA
2006-4-20-1R	R.H. SIDE WALL Bend Location: A--	1.00	EA
2006-4-30-3	INNER LINER R.H. SIDE 18 GA PTLOK	1.00	EA
2006-4-30-4	INNER LINER L.H. SIDE Bend Location: B--	1.00	EA
2006-4-300-5	CONDENSING UNIT PLATE 18 GA. GALV., 32 5/8 X 17 7/8"	1.00	EA
2006-4-40-1	FRAME END Bend Location: B--	4.00	EA
2006-4-40-2	FRAME SIDE Bend location: B--	4.00	EA
2006-4-60-10	HORIZONTAL GASKET SUPPORT 18 GA. PTLOK	2.00	EA
2006-4-60-10	HORIZONTAL GASKET SUPPORT 18 GA. PTLOK	2.00	EA
2006-4-60-11	VERTICAL GASKET SUPPORT 18 GA. PTLOK	2.00	EA
2006-4-60-11	VERTICAL GASKET SUPPORT 18 GA. PTLOK	2.00	EA
2018-8-600-2	LONG COPPER TUBE 1/2" OD COPPER, 18" LENGTH	2.00	EA
2019-18-100-2	MOUNTING BRACKET Bend Location: B--	1.00	EA
2020-2-100-5	36 SOLID-STATE RELAY MOUNTING PLATE Bend Location: B--	1.00	EA
2020-2-100A	HINGED INTELLUS BRACKET WELD	1.00	EA
2020-2-100B	CONTROL BOX LOCK CATCH PLATE WELD	1.00	EA
2020-2-80-1L	LH CONTROL TOP CORNER Bend Location: B--	1.00	EA
2020-2-80-1R	RH CONTROL TOP CORNER 18 GA PTLOK	1.00	EA
2020-3-100-1	41 CONTROL BOX DOOR W/ TABLET Bend Location: A--	1.00	EA
2020-3-100-2	41 CONTROL BOX ENCLOSURE Bend Location: A--	1.00	EA
2020-3-100-3	41 CONTROL BOX TOP COVER 18 GA PTLOK	1.00	EA
2020-3-100-4	41 CONTROL BOX DOOR STOP Bend Location: B--	1.00	EA
2020-3-80-2	41 CONTROL TOP REAR PANEL 22 GA PTLOK	1.00	EA
2021-8-170-1	HEATER BRACKET	2.00	EA
70-1-70-8A	HORIZ. LIGHT FIXTURE CLIP .050 ALUM. 3/4" x 1 1/4"	4.00	EA
70-2-300-16	SUCTION ACCUMULATOR MTG. BRACKET Bend Location: B--	2.00	EA
95-3-120-3	SIDE PANEL CHANNEL Bend Location: B--	2.00	EA
95-3-120-3	SIDE PANEL CHANNEL Bend Location: B--	2.00	EA
95-3-20-7A	GUSSET 18 GA. PTLOK, 8 X 8" (EACH), 28 13/32 X 40 3/8" (BATCH 32)	4.00	EA
95-3-20-8L	REAR CORNER LINER SUPPORT Bend Location: B--	1.00	EA
95-3-20-8R	REAR CORNER LINER SUPPORT Bend Location: B--	1.00	EA

**Partial Parts List For Job #: 21244-0000**

Date Printed: March 13, 2024

For Sales Order: Item: LT41VLSWC8T

Part #	Description	QTY	UOM
95-3-300-5	TOP LH PANEL <i>Bend Location: B--</i>	1.00	EA
95-3-300-6	TOP RH SCREEN <i>Bend Location: B--</i>	1.00	EA
95-3Q-500-1	SKID SHIM 2 X 4 LUMBER 15 1/2" LENGTH	2.00	EA
95-4-50-5	CORNER-DOOR CORE <i>Bend Location: FLAT</i>	4.00	EA
95481	SPORLAN SOLENOID VALVE MOUNT 16 GA GALV 1 5/8" X 3 3/8"	4.00	EA
CMS-044 (MACH.	6-32 X 1/2" PHILLIPS ROUND HD. S.S.	4.00	EA
CMS-069 (MACH.	10-24 UNC X 3/8, PH. RD. HD., S.S. 10-24 UNC X 3/8" PHILLIPS, ROUND HEAD, STAINLESS-STEEL MACHINE SCREW	2.00	EA
CMS-074 (MACH.	6-32 X 1" PHILLIPS RD. HD., S.S.	4.00	EA
CMS-111 (MACH.	6-32 X 1/4" PHIL. PAN HEAD, ZINC LED LIGHT MOUNTING SCREWS. Thrd-Cutting Screw for Metal&Plastic (Type F) Pan Head Phil, Zinc-Plated Steel, 6-32 Thread, 1/4" L	31.00	EA
CMS-111 (MACH.	6-32 X 1/4" PHIL. PAN HEAD, ZINC LED LIGHT MOUNTING SCREWS. Thrd-Cutting Screw for Metal&Plastic (Type F) Pan Head Phil, Zinc-Plated Steel, 6-32 Thread, 1/4" L	31.00	EA
CSM-004 (S.M. SCREW)	6 X 1/2" NICKEL PL. PHIL. TRUSSHEAD #6128-148N SERVICE SUPPLY CO.	4.00	EA
CSM-004 (S.M. SCREW)	6 X 1/2" NICKEL PL. PHIL. TRUSSHEAD #6128-148N SERVICE SUPPLY CO.	4.00	EA
CSM-039 (S.M. SCREW)	6 X 1/2" SELF DRILLING #410 S.S. HEX HEAD, CORROSION RESISTANT COATING	15.00	EA
CWN-002 (WELDNUT)	10-24 UNC DUAL SPUR SPOTWELD NUT BUCKEYE ND-1514; 10-24 UNC; DUAL SPUR; REQUIRES A Ø 1/4" PILOT HOLE; C1010 LOW-CARBON STEEL.	4.00	EA
CWN-002 (WELDNUT)	10-24 UNC DUAL SPUR SPOTWELD NUT BUCKEYE ND-1514; 10-24 UNC; DUAL SPUR; REQUIRES A Ø 1/4" PILOT HOLE; C1010 LOW-CARBON STEEL.	4.00	EA
CWN-002 (WELDNUT)	10-24 UNC DUAL SPUR SPOTWELD NUT BUCKEYE ND-1514; 10-24 UNC; DUAL SPUR; REQUIRES A Ø 1/4" PILOT HOLE; C1010 LOW-CARBON STEEL.	4.00	EA
CWN-003 (WELDNUT)	1/4-20 UNC, DUAL SPUR, C1010 STEEL BUCKEYE ND-2118, DUAL-SPUR, 1/4-20 UNC SPOTWELD NUT; C1010 LOW-CARBON STEEL, REQUIRES A Ø 5/16" PILOT HOLE.	8.00	EA
CWN-003 (WELDNUT)	1/4-20 UNC, DUAL SPUR, C1010 STEEL BUCKEYE ND-2118, DUAL-SPUR, 1/4-20 UNC SPOTWELD NUT; C1010 LOW-CARBON STEEL, REQUIRES A Ø 5/16" PILOT HOLE.	8.00	EA
CWN-007	6-32 STEEL WELD NUT 0.185" HOLE SIZE, 0.18" TO 0.19" OD BARREL SIZE.	2.00	EA
CWN-007	6-32 STEEL WELD NUT 0.185" HOLE SIZE, 0.18" TO 0.19" OD BARREL SIZE.	4.00	EA
CWS-023 (DRYWALL ...	8 X 3 DRYWALL SCREW BUGLE HD DRYWALL SCREW ZINC/YELLOW PHILIPS COARSE THREAD	20.00	EA
CWS-024 (DRYWALL ...	8 X 2 PARTIAL BOARD SCREW PARTIAL THREAD 8 X 2 PHILIPS HEAD FLAT HEAD WITH NIBS, BLACK PHOS STEEL	46.00	EA
EBR-146 (BREAKER)	15A, 1-P, ALTECH 1DU15L, UL489 ALTECH CORP. CATALOG #1DU15L, DIN MOUNT (AC) CIRCUIT BREAKER, 15A, 1-POLE, 10kA @ 120/240/277VAC, 50/60Hz, L-SERIES; SCREW TERMINALS, GREEN: OFF/TRIP, RED: ON, LINE/LOAD REVERSIBLE, HIGH MAGNETIC TRIP POINT, UL 489 & CE	1.00	EA
EBR-147 (BREAKER)	2A, 1-P, ALTECH 1DU2L, UL489 ALTECH CORP. CATALOG #1DU2L, DIN MOUNT (AC) CIRCUIT BREAKER, 2A, 1-POLE, 10kA @ 120/240/277VAC, 50/60Hz, L-SERIES; SCREW TERMINALS, GREEN: OFF/TRIP, RED: ON, LINE/LOAD REVERSIBLE, HIGH MAGNETIC TRIP POINT, UL 489 & CE	3.00	EA
EBR-163 (BREAKER)	3A, 1-P, ALTECH 1DU3L, UL489 ALTECH 1DU3L, DIN MOUNT (AC) CIRCUIT BREAKER, 3A, 1-POLE, 10kA @ 120/240/277VAC, 50/60Hz, L-SERIES; SCREW TERMINALS, GREEN: OFF/TRIP, RED: ON, LINE/LOAD REVERSIBLE, HIGH MAGNETIC TRIP POINT, UL 489 & CE	2.00	EA
EBR-180 (BREAKER)	25A, 1-P, ALTECH 1DU25L, UL489 ALTECH 1DU25L, DIN MOUNT (AC) CIRCUIT BREAKER, 25A, 1-POLE, 10kA @ 120/240/277VAC, 50/60Hz, L-SERIES; SCREW TERMINALS, GREEN: OFF/TRIP, RED: ON, LINE/LOAD REVERSIBLE, HIGH MAGNETIC TRIP POINT, UL 489 & CE	1.00	EA
EBU-016 (BUSHING)	SR7N-2 HEYCO INSULATOR BUSHING HEYCO #1244 STRAIN RELIEF BUSHING; BLACK, UL/CSA STRAIGHT THRU-TYPE; FOR USE W/ Ø 0.430" MAX. ROUND CABLES; FOR USE W/ (EPL-030) 7-WIRE CORDSET. REQUIRES A 0.66" WIDE X 0.75" LONG OVAL MOUNTING HOLE.	1.00	EA
EBU-016 (BUSHING)	SR7N-2 HEYCO INSULATOR BUSHING HEYCO #1244 STRAIN RELIEF BUSHING; BLACK, UL/CSA STRAIGHT THRU-TYPE; FOR USE W/ Ø 0.430" MAX. ROUND CABLES; FOR USE W/ (EPL-030) 7-WIRE CORDSET. REQUIRES A 0.66" WIDE X 0.75" LONG OVAL MOUNTING HOLE.	1.00	EA
EBU-016 (BUSHING)	SR7N-2 HEYCO INSULATOR BUSHING HEYCO #1244 STRAIN RELIEF BUSHING; BLACK, UL/CSA STRAIGHT THRU-TYPE; FOR USE W/ Ø 0.430" MAX. ROUND CABLES; FOR USE W/ (EPL-030) 7-WIRE CORDSET. REQUIRES A 0.66" WIDE X 0.75" LONG OVAL MOUNTING HOLE.	1.00	EA
EBU-016 (BUSHING)	SR7N-2 HEYCO INSULATOR BUSHING HEYCO #1244 STRAIN RELIEF BUSHING; BLACK, UL/CSA STRAIGHT THRU-TYPE; FOR USE W/ Ø 0.430" MAX. ROUND CABLES; FOR USE W/ (EPL-030) 7-WIRE CORDSET. REQUIRES A 0.66" WIDE X 0.75" LONG OVAL MOUNTING HOLE.	1.00	EA
ECL-004 (CLAMP)	3/16" NYLON CABLE CLAMP WHITE	2.00	EA
ECO-003 (CONNECTOR)	T&B 1/2" STRAIGHT, FLEXIBLE STCXC270, 1/2" SQUEEZE FLEX CONNECTOR; UL- & CSA-APPROVED.	1.00	EA
ECO-061 (CONNECTOR)	WAGO 231-304/026-000 FOUR CIRCUIT WAGO #231-304/026-000; 4-POLE ORANGE FEMALE CONNECTOR W/ CAGE CLAMP, FOR USE W/ AWG #28-12 WIRE, 300 V/15 A MAX., UL & CSA APPROVED.	2.00	EA
ECO-121 (CONNECTOR)	WAGO 5 CONDUCTORS LEVER CONNECTOR WAGO 221-415, 24-12 GA. WIRE, 600V, 20A, 3 COMMON TERMINALS, UL, CE, EN, 25/BOX	5.00	EA
EFA-039 (FAN)	4YD79 AXIAL FAN CORD SET GRAINGER #4YD79; CORD SET REQUIRED FOR USE OF GRAINGER #5C155 (EFA-038) AXIAL FAN; UL-APPROVED (#E88231) & CSA-APPROVED (#LL54496).	2.00	EA
EFA-039 (FAN)	4YD79 AXIAL FAN CORD SET GRAINGER #4YD79; CORD SET REQUIRED FOR USE OF GRAINGER #5C155 (EFA-038) AXIAL FAN; UL-APPROVED (#E88231) & CSA-APPROVED (#LL54496).	2.00	EA
EFA-045 (FAN)	4YD87 AXIAL FAN GUARD	2.00	EA

**Partial Parts List For Job #: 21244-0000**

Date Printed: March 13, 2024

For Sales Order: Item: LT41VLSWC8T

Part #	Description	QTY	UOM
EFA-045 (FAN)	GRAINGER 4YD87; FITS AXIAL FANS W/ 4-1/8" X 4-1/8" O.C. MOUNTING HOLE DIMENSIONS. 4YD87 AXIAL FAN GUARD	2.00	EA
EFA-091 (FAN)	GRAINGER 4YD87; FITS AXIAL FANS W/ 4-1/8" X 4-1/8" O.C. MOUNTING HOLE DIMENSIONS. AXIAL FAN, 230VAC, 31CFM, 50/60HZ	1.00	EA
EFA-111 (FAN)	ORION OA80AP-22-1WB, SHADED-POLE MOTOR, 0.06a @ 230VAC, 12W, 31 cfm @ 50/60Hz, 3000 RPM, 29 dB, -20°-80°C TEMP. OPERATING RANGE, CW ROTATION, SUPPLIED W/ TWO (2) 22 GA. X 12" LEAD WIRES, DUAL BALL BEARING, RoHS-COMPLIANT, 3.15" X 3.15" X 1.50"	2.00	EA
EFA-111 (FAN)	OA125AP-22-1TB ORION FAN 220V 75CFM	2.00	EA
EFA-111 (FAN)	OA125AP-22-1TB AC, 220V; 120x120x25mm; Sq; 75CFM; 14W; 42dBA; 2900RPM; Terminals	2.00	EA
EHE-046 (HEATER)	250W @ 230V U-TUBE HEATER, RED OPEN AIR HEATER, 250W @ 230V, MARKED WITH RED PAINT, WATLOW PART RBN36B1H-2, 40" OAL X 0.315" DIA., FORMED IN FIGURE 1, A = 17.5625"; R = 1"; W/ EPOXY SEAL AND 36" SILICONE WIRE LEADS W/ SHRINK BOOT EPOXY SEAL, QUOTE # 710091 See drawing number 2006-1-300-7	1.00	EA
EHE-052 (HEATER)	230V, 11W, FLEX. BRAID ALUM HEATER 11W, 230VAC FLEXIBLE BRAID ALUMINUM HEATER, 8" LONG, REFRIGERATION HARDWARE SUPPLY CORP NO. 64-402	2.00	EA
EHE-068 (HEATER)	230V, 51W, FLEX. BRAID ALUM HEATER 51W, 230VAC, 0.22A, FLEXIBLE BRAID ALUMINUM HEATER, 40" LONG, REFRIGERATION HARDWARE SUPPLY CORP NO. 64-416	2.00	EA
ELD-033	DALI TO PWM CONVERTER (ONE DALI ) MEAN WELL DAP-04	1.00	EA
ELD-066	12 X 24 WHITE LED TILE W/COATING 12" x 24" LED tile with conformal coating. With Samsung LM301B Mid Power LED, 4000K with CR190 SH flux bin. 12 STRINGS OF LEDs WITH 13 LEDs PER STRING (156).	3.00	EA
ELD-066	12 X 24 WHITE LED TILE W/COATING 12" x 24" LED tile with conformal coating. With Samsung LM301B Mid Power LED, 4000K with CR190 SH flux bin. 12 STRINGS OF LEDs WITH 13 LEDs PER STRING (156).	3.00	EA
ELD-071	LED DRIVER MEANWELL HLG-185H-C1400B MEANWELL HLG-185H-C1400B LED DRIVER, CONSTANT CURRENT, 1400mA OUTPUT CURRENT, 200.2W, 71-143V OUTPUT, 90-305VAC INPUT VOLTAGE, WITH DIMMING, -40 TO 90°C OPERATING TEMP RANGE, RoHS, 8.03"L x 2.68"W x 1.53"H, IP65 RATING, 3 IN 1 DIMMING	1.00	EA
ELD-071	LED DRIVER MEANWELL HLG-185H-C1400B MEANWELL HLG-185H-C1400B LED DRIVER, CONSTANT CURRENT, 1400mA OUTPUT CURRENT, 200.2W, 71-143V OUTPUT, 90-305VAC INPUT VOLTAGE, WITH DIMMING, -40 TO 90°C OPERATING TEMP RANGE, RoHS, 8.03"L x 2.68"W x 1.53"H, IP65 RATING, 3 IN 1 DIMMING	1.00	EA
EPC-032 (USB CABLE)	PANEL MOUNT USB CABLE, A-A F/M, 3FT STARTECH USBPNLAFAM3	1.00	EA
EPC-044 (COMPUTER)	SIBO, 10" ANDROID TOUCHSCREEN SHENZHEN SIBO INDUSTRIAL & DEVELOPMENT CO, LTD. MODEL Q8919 TOUCHSCREEN HMI WITH:	1.00	EA
EPC-057 (ETHERNET)	ETHERNET CABLE - 5' CAT5 CABLE BLK 10X6-12205, 5 FT., BLACK, UL-APPROVED, RoHS-COMPLIANT, WITHOUT PVC BOOT	1.00	EA
EPL-004 (PLUG)	3 TERM., MOLEX 1396, 03-09-2031 Conn Housing PL 3 POS 5.03mm Crimp ST Panel Mount Bag, W/MTG. EARS	4.00	EA
EPL-006 (PLUG)	MALE PLUG W/LEAD 18/3 W/54"LD WHITE EMF CORP. M-3-7/PM0030; 3-WIRE CORDSET, W/ 54" LENGTH AWG #18 STRANDED Cu WIRES & 300 V/105°C MAX., 1/2" PRESTRIPPED WHITE PVC INSULATION; UL- & CSA-RECOGNIZED.	1.00	EA
EPL-041 (PLUG)	9 TERM., AMP MATE-N-LOK 1-480706-0 TYCO UNIVERSAL MATE-N-LOK 1-480706-0; 9-TERMINAL PLUG, 19 A @ 600 V MAX., 105°C MAX.; USE W/ (ERE-046); UL- & CSA-APPROVED.	3.00	EA
EPL-052 (PLUG)	EMF CORP. P-7-1 CORDSET (7-COLOR) EMF CORP. P-7-1 7-WIRE CORDSET; SEVEN COLORS, SEVEN (7) AWG 18 X 48" STRANDED Cu WIRES, SEVEN-PIN RIGHT-ANGLE PLUG MOLDED ON 3/8" X 44" LENGTH WHITE PVC INSULATION, 4" LENGTH EXPOSED 1/4" PRESTRIPPED WIRE LEADS.	1.00	EA
EPL-052 (PLUG)	EMF CORP. P-7-1 CORDSET (7-COLOR) EMF CORP. P-7-1 7-WIRE CORDSET; SEVEN COLORS, SEVEN (7) AWG 18 X 48" STRANDED Cu WIRES, SEVEN-PIN RIGHT-ANGLE PLUG MOLDED ON 3/8" X 44" LENGTH WHITE PVC INSULATION, 4" LENGTH EXPOSED 1/4" PRESTRIPPED WIRE LEADS.	1.00	EA
EPL-059 (PLUG)	TYCO 172167-1 MINI M-N-L , 4 SOCKET 172167-1 MINI MATE-N-LOK Connector, Soft Shell; 4; 600 VAC; Plug; Double-Row; Nylon; Plug; Natural	1.00	EA
EPS-031 (POWER SU...)	120/240V IN, 12VDC OUT, 75W MEAN WELL LRS-75-12	1.00	EA
ERE-012 (RECEPTACLE)	3 HOLE W/EARS AND TABS PLUG & SOCKET HOUSING, RCPT, NYLON	4.00	EA
ERE-014 (RECEPTACLE)	FEMALE PLUG W/LEAD, 18/3W/12"LD EMF CORP. PM0020; 3-WIRE FEMALE RECEPTACLE, W/ THREE (3) 12" LENGTH X AWG #18 STRANDED Cu WIRES & 300 V/105°C MAX., 1/2" PRESTRIPPED WHITE PVC INSULATION; REQUIRES A 1" X 1" SQUARE MOUNTING HOLE; FITS A 0.060" MAX. PANEL THICKNESS; UL- & CSA-RECOGNIZED.	1.00	EA
ERE-083 (RECEPTACLE)	R-7-1, 108" SOLID WIRE 7-COLOR SEVEN SOCKET / SEVEN COLOR, 108" SOLID WIRE LEN., W / WIRE STRIPPED AT ENDS	2.00	EA
ERL-030 (RELAY)	DPDT, 30A., 240VAC	1.00	EA
ERL-060 (RELAY)	S.S., N.O., 3-32 VDC INPUT, 5A, DIN CARLO GAVAZZI RP1A23D5M1; N.O., 3-32 VDC INPUT, 24-240 VAC OUTPUT, 5 A MAX. LOAD CURRENT, ZERO-SWITCH TRIAC; DIN MOUNT, PUSH-TO-RELEASE WIRING CONNECTIONS; UL- & CSA-APPROVED.	7.00	EA
ERL-063 (RELAY)	P&B/TE, DPDT, 30A., 240VAC POTTER&BRUMFIELD/TE CONNECTIVITY, T92S11A22-240 Power Relay, DPDT, 240 VAC, 30 A, T92 Series, Panel Mount, Non Latching	4.00	EA
ESW-130 (SWITCH)	DPST ROCKER SWITCH, 20A, 125V PRCDA1-20F-BB0BW DPST ROCKER SWITCH, ON-OFF, 20A @ 125V, 10A @ 250V, PANEL MOUNT, SNAP-IN, BLACK, I-O MARKINGS, RoHS COMPLIANT,	1.00	EA
ETC-052 (TEMP. CON...)	3 WIRE RTD, 100 OHMS A/100-3W-2.5"-3/16"-12"CL2P-HS, 144" LEADS, 3 WIRE RTD SENSOR, 100 OHMS	2.00	EA
ETC-108 (TEMP. CON...)	INTELLUS CONTROL OVERLAY - CHINA INTELLUS CONTROL OVERLAY WITH GREEN BACKGROUND FOR CHINA	1.00	EA
ETC-109 (TEMP. CON...)	INTELLUSULTRA GEN2, LOW END VER2 Watlow PSIB-ASIC-0005	1.00	EA
ETE-017 (TERMINAL)	WAGO 2004-1401, GRAY, 4 TERMINAL 5 mm TERMINAL BLOCK WIDTH FOR MOUNTING ON DIN RAIL	12.00	EA
ETE-018 (TERMINAL)	WAGO 2004-1404, BLUE, 4 TERMINAL 5 mm TERMINAL BLOCK WIDTH FOR MOUNTING ON DIN RAIL	6.00	EA

Partial Parts List For Job #: 21244-0000

Date Printed: March 13, 2024

For Sales Order: Item: LT41VLSWC8T

Part #	Description	QTY	UOM
ETE-019 (TERMINAL)	WAGO 2004-1407, GROUND, 4 TERMINAL 5 mm TERMINAL BLOCK WIDTH FOR MOUNTING ON DIN RAIL	3.00	EA
ETE-020 (TERMINAL)	WAGO 2004-1491, END PLATE 5 mm TERMINAL BLOCK WIDTH FOR MOUNTING ON DIN RAIL	3.00	EA
ETE-021 (TERMINAL)	WAGO 2004-402; GRAY TERM. JUMPER WAGO #2004-402; LIGHT GRAY, 2-PIN PUSH-IN JUMPER BAR; FOR USE W/ 5 mm WIDE TERMINAL BLOCKS; 32 A MAX.	7.00	EA
ETE-022 (TERMINAL)	WAGO 249-116, END TERMINAL ANCHOR FOR 5 mm TERMINAL BLOCK WIDTH	8.00	EA
ETE-030 (TERMINAL)	16-CONDUCTOR, 4-TERMINAL SNAP-IN WAGO #862-0594, BLACK, NO LABELS. MTG REQUIRES TWO Ø5.2mm HOLES SPACED 12mm, PACKS OF 200; UL- & CUL-APPROVED.	1.00	EA
ETE-032 (TERMINAL ...)	JUMPER FOR WAGO TERM BLOCK ETE-030 WAGO #281-482, GRAY, 2-PIN PUSH-IN JUMPER; UL- & CSA-APPROVED. 100 PC. PER BOX	1.00	EA
GDP-012 (GLASS DBL...)	230V HEATED DBL PANE GLASS W/ARGON ARGON FILLED, 49 3/4 X 22 3/16 X 3/4", TEMPERED, DOUBLE-PANE, 230VAC, 6.41 Watts/Ft <sup>2</sup> WITH DIODE, SUPPLIED BY CIG GLASS ( <a href="http://www.cigglass.com/">http://www.cigglass.com/</a> ). ORDER WITH APPROX. 24" LEAD WIRES EXITING AT MIDDLE OF SHORT EDGE OF GLASS; DIODE TO BE EXTERNALLY ATTACHED TO LEAD WIRE FOR EASY REMOVAL IF NECESSARY.	2.00	EA
HBS-001 (SPACER)	FEMALE THREADED CIRCUIT BOARD POST MICRO PLASTICS INC #27FTP00250	6.00	EA
HDR-007 (DRAIN)	328 005 THRU HULL DRAIN 3/4" W/NUT	1.00	EA
HHA-008 (HANDLE)	BLACK SOUTHCO P1-30-101-11 ABS LOW GLOSS	1.00	EA
HHI-029 (HINGE)	SOUTHCO EH-5A-4G4-50 SERIES 40 BLACK POWDER COATED MOUNTING METHOD: M5 FOUR THREADSTUDS AT .984" BETWEEN CENTER.	2.00	EA
HLL-001 (LEVELER)	1/2"-2.264" OVERALL HEIGHT SK-9246-FT 1/2-13 X 2 1/4 LG LVL R	4.00	EA
HLO-019 (LOCK)	COMPRESSION LOCK (1/4 TURN) 040-1-1-1-42-59, ELECTRAUSA chrome with slotted lock key	1.00	EA
HMA-010	S-MAGNET 1" X 1/2" X 1/8" THICK BC82CS-S, 1" X 1/2" X 1/8" THICK WITH COUNTERSUNK HOLE FOR #6 SCREW, RECTANGULAR BLOCK MAGNET, N42, NICKEL PLATING, MAGNETIZED WITH THE SOUTH POLE ON THE COUNTERSUNK SIDE	1.00	EA
HNU-052 (NUT)	M5, .8 COARSE THREAD NUT	8.00	EA
HNU-076 (RIVET NUT)	6-32 RUBBER-INSULATED RIVET NUT McMaster-Carr 93495A120, Rubber-Insulated Brass Rivet Nut, 6-32 Thread, 1/64" - 5/32" Material Thick Range	4.00	EA
HST-002 (PORT PLUG)	1 1/4" ACCESS PORT PLUG SHAW AGRO DEVELOPMENT #51013	1.00	EA
LNA-015 (NAILS)	S12D-FH 3 1/4" FRAMING NAILS BOSTITCH S12D-FH Framing Nails, 28 Degree, Wire Weld, 3-1/4-Inch x .120-Inch, 2000-Pack	15.00	EA
LST-006 (STAPLES)	1/2 CROWN X 1 3/4 LEG 16GA. STAPL INTERCHANGE GS16-134 GALVANIZED SHEATHING STAPLE	22.00	EA
LST-007 (STAPLES)	15/16 CROWN X 3/4 LEG 16GA. STAPL ROOFING STAPLE	144.00	EA
LST-008 (STAPLES)	15/16 CROWN X 1 LEG 16GA. STAPLE INTERCHANGE WC GALVANIZED STAPLE	287.00	EA
RAC-003 (SUCTION A...)	3702 REFRIG. RESEARCH, 5/8 ODF REFRIGERATION RESEARCH #3702, 5/8" ODF SOLDER CONNECTIONS, 4.1 lb R-22 HOLDING CAPACITY, 4.2 lb R-134a HOLDING CAPACITY	1.00	EA
RBR-001 (RETURN BEND)	1/2" O.D. X .028" W X 1-3/8" #27826 WROT COPPER, MALE ENDS	2.00	EA
RBR-001 (RETURN BEND)	1/2" O.D. X .028" W X 1-3/8" #27826 WROT COPPER, MALE ENDS	2.00	EA
RCO-084 (UNIT COOLER)	BOHN TA23SG W/ 230/1/60 FANS BOHN, 230/1/60 Hz, 1.2 AMPS w/ 3 fan motors	2.00	EA
RCU-375 (CONDENSIN...)	M6KP-0102-CAZ-372.449A 220-240/1/50 EMERSON M6KP-0102-CAZ-372; 1 HP, AIR-COOLED, R-449A, 220/240-1-50, 3/8" LIQUID-LINE SWEAT CONNECTION, 1/2" SWEAT SUCTION-LINE CONNECTION, 1/4" DISCHARGE LINE, COMP. RLA: 9.4A, FAN RLA: 0.76A, TOT RLA: 10.16A	1.00	EA
RES-002 (ELBOW)	ELBOW, 1/2"NOM, 5/8"OD WROT, 90° STREET FITTING W-02326	1.00	EA
RFN-001 (FLARE NUT)	FREY 14FS4, 1/4" SHORT TYPE FORGED	2.00	EA
RFN-003 (FLARE NUT)	NS4-8 1/2" SAE FLARE	2.00	EA
RPR-006 (PRES. REG...)	SPORLAN CROT-6-30/110, 5/8" ODF SPORLAN CROT-6-30/110 5/8" ODF; 30-110 psig PRESSURE RANGE; 5/8" ODF SOLDER CONNECTIONS; SUPPLIED W/ INLET STRAINER.	1.00	EA
RTE-019 (TEE)	TEE, WROT 1/4"NOM SIZE, W-04000 (3/8" OD COPPER TUBE)	6.00	EA
RTE-021 (TEE)	TEE, WROT 1/2"OD, 3/8" NOM., W-04001 MUELLER BRASS W-04001	1.00	EA
RTE-030 (TEE)	W-04002, 1/2" X 1/2" X 3/8"	4.00	EA
RVH-002 (HAND VALVE)	EBVP-1030 3/8" ODF BALL VALVE WITH ACCESS PORT, PARKER 502054	1.00	EA
RVQ-002 (BODY)	Q BODY 3 X 4 SAE 168000, 3/8 SAE INLET, 1/2 SAE OUTLET	2.00	EA
RVQ-101 (CARTRIDGE)	QC-1 YELLOW Q CARTRIDGE/ID TAG 168072	2.00	EA
RVQ-201 (ELEMENT)	KT-43-VC 30' ELEMENT 180269	2.00	EA
RVS-027 (VALVE)	FEMALE DIN CONNECTOR ASSY GOES W/RVS-028 ON RVS-030 & RVS-031	4.00	EA
RVS-028 (VALVE)	OMKC-1E 220-240/50 COIL	4.00	EA

Partial Parts List For Job #: 21244-0000

Date Printed: March 13, 2024

For Sales Order: Item: LT41VLSWC8T

Part #	Description	QTY	UOM
RVS-031 (VALVE)	SOLENOID VALVE XUP-7, LESS COIL SPORLAN XUP, W/ EXTENDED 3/8" ODF SOLDER CONNECTIONS	4.00	EA
SCL-005 (CLIP)	#66 STAINLESS STL CLIP 0066-000008	8.00	EA
SPI-004 (PILASTER)	48" 0060-007048 60 STANDARD	4.00	EA
SWS-115 (SHELF)	SHELF, 36 1/4" X 27" W/WHITE EPOXY-POLYESTER HYBRID FINISH, DWG. SWS-115(SOLIDWORKS), 41 SERIES	2.00	EA
XGA-051 (GASKET)	39" X 59 3/4" OD, 4 SIDED BLACK 02-126 BLACK MAGNETIC DOOR GASKET, 4 SIDED	1.00	EA
XGA-053 (GASKET)	54 5/8" X 28 1/8" O.D., 02-126 4 SIDED, BLACK	1.00	EA
XGA-053 (GASKET)	54 5/8" X 28 1/8" O.D., 02-126 4 SIDED, BLACK	1.00	EA
XGR-002 (GROMMET)	3/4" I.D. X 7/8" HOLE	12.00	EA
XGR-002 (GROMMET)	3/4" I.D. X 7/8" HOLE	2.00	EA
XGR-003 (GROMMET)	1-1/4 ID X 1-1/2 HOLE LAVELLE #1047 LAVELLE #1047; 1-7/8" O.D. X 1-1/4" I.D. X Ø 1-1/2" PANEL HOLE X 3/8" O.A. THICKNESS X 1/16" GROOVE WIDTH, BLACK RUBBER.	1.00	EA
XGR-003 (GROMMET)	1-1/4 ID X 1-1/2 HOLE LAVELLE #1047 LAVELLE #1047; 1-7/8" O.D. X 1-1/4" I.D. X Ø 1-1/2" PANEL HOLE X 3/8" O.A. THICKNESS X 1/16" GROOVE WIDTH, BLACK RUBBER.	1.00	EA
XGR-003 (GROMMET)	1-1/4 ID X 1-1/2 HOLE LAVELLE #1047 LAVELLE #1047; 1-7/8" O.D. X 1-1/4" I.D. X Ø 1-1/2" PANEL HOLE X 3/8" O.A. THICKNESS X 1/16" GROOVE WIDTH, BLACK RUBBER.	2.00	EA
XGR-005 (GROMMET)	3/8" I.D. X 1/2" HOLE, BLACK VINYL BLACK VINYL, G1089 PLASTIC	2.00	EA
XGR-011 (GROMMET)	1/8ID X 11/32OD X 7/32T X 1/4 HOLE McMASTER-CARR #9600K11; 1/8" I.D. X 11/32" O.D. X 7/32" THICK X 1/16" GROOVE WIDTH, REQUIRES A Ø 1/4" MOUNTING HOLE; SBR (STYRENE BUTADIENE RUBBER) MATL; 50 SHORE A DUROMETER HARDNESS, -35° F TO + 180° F TEMP. RANGE.	2.00	EA
XGR-022 (GROMMET)	WHT LCKING GROMMET FOR 3/8 DIA HOLE McMaster-Carr 5269T232, Locking Grommet for 3/8" Hole Diameter, 3/16" ID, White	3.00	EA
XGR-022 (GROMMET)	WHT LCKING GROMMET FOR 3/8 DIA HOLE McMaster-Carr 5269T232, Locking Grommet for 3/8" Hole Diameter, 3/16" ID, White	3.00	EA
YDR-002	ACCESS PORT SLEEVE, 1 1/4" ID COMPONENT HARDWARE R90-1002 Flexible Insulating Conduit Sleeve	1.00	EA
YHI-002 (HINGE)	IN LINE KNUCKLE STYLE HINGE SOUTHCO 96-10-400-11	5.00	EA
YHI-002 (HINGE)	IN LINE KNUCKLE STYLE HINGE SOUTHCO 96-10-400-11	4.00	EA
YHI-002 (HINGE)	IN LINE KNUCKLE STYLE HINGE SOUTHCO 96-10-400-11	4.00	EA
YHP-002 (HOLE PLUG)	3/4" HOLE WHITE DP-750 PLASTIC	1.00	EA
YHP-012 (HOLE PLUG)	1/2" WHITE DP-500 PLASTIC 2644	2.00	EA
YHS-001 (SPACER)	HINGE SPACER 95-3-50-5A 30% GLASS FILLED NYLON/ BLACK	5.00	EA
YHS-001 (SPACER)	HINGE SPACER 95-3-50-5A 30% GLASS FILLED NYLON/ BLACK	4.00	EA
YHS-001 (SPACER)	HINGE SPACER 95-3-50-5A 30% GLASS FILLED NYLON/ BLACK	4.00	EA
YWA-001 (WASHER)	.040 PLASTIC SHOULDER WASHER 24SW0006 0.040" Nylon 6/6 Shoulder Washer with Natural Finish, Off White	24.00	EA
YWA-001 (WASHER)	.040 PLASTIC SHOULDER WASHER 24SW0006 0.040" Nylon 6/6 Shoulder Washer with Natural Finish, Off White	24.00	EA