

Partial Parts List For Job #: 21402-0000

Date Printed: March 21, 2024

For Sales Order: 33551 Item: E41L2C8TSWH3

Part #	Description	QTY	UOM
2004-1-300-27	SINGLE TERMINAL BLOCK MTG BRACKET 18 GA. GALV. 3 9/16" X 1 29/32"	1.00	EA
2006-1-300-6	HEATER BRACKET 22 GA GALV	4.00	EA
2006-1-50-1	DOOR FACE PANEL Bend Location: A--	1.00	EA
2006-1-50-2	DOOR CORE Bend Location: A--	1.00	EA
2006-1-50-3	CU-41 DOOR INSULATION POLYSTYRENE, 1 7/8 X 39 X 59 1/2"	1.00	EA
2006-5-150-1	REAR WALL PLENUM 22 GA PTLOK	1.00	EA
2006-5-150A	REAR WALL PLENUM SLIDER WELD ASS'Y	11.00	EA
2006-5-300-5	DRIP DEFLECTOR .040 ALUM., 30 3/4 x 1 1/2	1.00	EA
2010-1-120-1	41 LAMP BANK BASE 22 GA PTLOK	2.00	EA
2010-1-120-2	41 LAMP BANK COVER 22 GA PTLOK	2.00	EA
2015-2C-30	41 SERIES INSULATED BOX (PURCHASED) REV. B IS THE NEW BOX DESIGN OFF THE SALVAGNINI WITH NEW FLOOR DRAIN & BREAKER REDESIGN.	1.00	EA
2017-5-230-1	SHORT (NORMAL) SENSOR BRACKET Bend Location: B--	1.00	EA
2019-18-100-2	MOUNTING BRACKET Bend Location: B--	1.00	EA
2020-2-100-5	36 SOLID-STATE RELAY MOUNTING PLATE Bend Location: B--	1.00	EA
2020-2-100A	HINGED INTELLUS BRACKET WELD	1.00	EA
2020-2-100B	CONTROL BOX LOCK CATCH PLATE WELD	1.00	EA
2020-2-80-1L	LH CONTROL TOP CORNER Bend Location: B--	1.00	EA
2020-2-80-1R	RH CONTROL TOP CORNER 18 GA PTLOK	1.00	EA
2020-21-110-1	DRIVER RACK Bend Location: B--	1.00	EA
2020-3-100-1	41 CONTROL BOX DOOR W/ TABLET Bend Location: A--	1.00	EA
2020-3-100-2	41 CONTROL BOX ENCLOSURE Bend Location: A--	1.00	EA
2020-3-100-3	41 CONTROL BOX TOP COVER 18 GA PTLOK	1.00	EA
2020-3-100-4	41 CONTROL BOX DOOR STOP Bend Location: B--	1.00	EA
2020-3-80-2	41 CONTROL TOP REAR PANEL 22 GA PTLOK	1.00	EA
2020-3-80-4	41 EXTENDED(2") MECH SECTION COVER 22 GA PTLOK	1.00	EA
2021-8-170-1	HEATER BRACKET 22 GA GALV	2.00	EA
57458	DRIER CLAMP .040 ALUM. 2 X 9 17/32	1.00	EA
70-1S-30-1	RH END PLATE .50 ALUM. 3 1/4 x 5 23/32	1.00	EA
70-1S-30-2	RH END PLATE .50 ALUM. 3 1/4 x 5 7/32	1.00	EA
70-1Y-30-1	WATER VALVE MOUNTING BRACKET 22 GA PTLOK 2" X 3 3/4"	1.00	EA
70-2-300-16	SUCTION ACCUMULATOR MTG. BRACKET 16 GA. GALV. 4" X 7"	2.00	EA
95-3-120-3	SIDE PANEL CHANNEL Bend Location: B--	2.00	EA
95-3-120-3	SIDE PANEL CHANNEL Bend Location: B--	2.00	EA
95-3-300-2	CONDENSING UNIT PLATE 18 GA GALV. 17 3/8 x 32 5/8	1.00	EA
95-3-300-5	TOP LH PANEL Bend Location: B--	1.00	EA
95-3-300-6	TOP RH SCREEN Bend Location: B--	1.00	EA
95-4-300-97	RH RECEPTACLE PANEL Bend location: A--	1.00	EA
95-4-50-5	CORNER-DOOR CORE Bend Location: FLAT	4.00	EA
95-4-600-1	DRAIN LINE COPPER TUBING (LOWER) 1/2" OD COPPER, .25" LENGTH	1.00	EA
95-4-600-2	DRAIN LINE COPPER TUBING (UPPER) 1/2" OD COPPER, .13" LENGTH	1.00	EA
95-5-30-3	HUMIDIFIER FLOAT SWITCH BRACKET 22 GA. #304-4 S.S.	1.00	EA
95-5-300-9 (EUR)	WATER VALVE MOUNTING BRACKET 16 GA. GALV.	1.00	EA

Partial Parts List For Job #: 21402-0000

Date Printed: March 21, 2024

For Sales Order: Item: E41L2C8TSWH3

Part #	Description	QTY	UOM
95-5-40-10	WATER VALVE MOUNTING BOX 22 GA PTLOK	1.00	EA
95-5-40-11	WATER VALVE MOUNTING BOX COVER 22 GA PTLOK	1.00	EA
95-5-700-2	VERTICAL WATER LINE CUT TO 5 1/4" LENGTH, USE YTU-077	1.00	EA
95481	SPORLAN SOLENOID VALVE MOUNT 16 GA GALV. 1 5/8 X 3 1/4	3.00	EA
98-2-300-7	LH RECEPTACLE PANEL 22 GA PTLOK	1.00	EA
CMS-007 (MACH.	8-32 X 3/8", RD.HEAD, PHILLIPS,ZINC	2.00	EA
CMS-027 (MACH.	1/4-20 UNC X 1/2", ROUND HD., ZINC ZINC-CHROMATE PLATED	10.00	EA
CMS-044 (MACH.	6-32 X 1/2" PHILLIPS ROUND HD. S.S.	4.00	EA
CMS-069 (MACH.	10-24 UNC X 3/8, PH. RD. HD., S.S. 10-24 UNC X 3/8" PHILLIPS, ROUND HEAD, STAINLESS-STEEL MACHINE SCREW	2.00	EA
CMS-069 (MACH.	10-24 UNC X 3/8, PH. RD. HD., S.S. 10-24 UNC X 3/8" PHILLIPS, ROUND HEAD, STAINLESS-STEEL MACHINE SCREW	3.00	EA
CMS-069 (MACH.	10-24 UNC X 3/8, PH. RD. HD., S.S. 10-24 UNC X 3/8" PHILLIPS, ROUND HEAD, STAINLESS-STEEL MACHINE SCREW	2.00	EA
CMS-069 (MACH.	10-24 UNC X 3/8, PH. RD. HD., S.S. 10-24 UNC X 3/8" PHILLIPS, ROUND HEAD, STAINLESS-STEEL MACHINE SCREW	22.00	EA
CMS-074 (MACH.	6-32 X 1" PHILLIPS RD. HD., S.S.	4.00	EA
CMS-111 (MACH.	6-32 X 1/4" PHIL. PAN HEAD, ZINC LED LIGHT MOUNTING SCREWS. Thrd-Cutting Screw for Metal&Plastic (Type F) Pan Head Phil, Zinc-Plated Steel, 6-32 Thread, 1/4" L	48.00	EA
CSM-004 (S.M. SCREW)	6 X 1/2" NICKEL PL, PHIL. TRUSSHEAD #6128-148N SERVICE SUPPLY CO.	4.00	EA
CSM-004 (S.M. SCREW)	6 X 1/2" NICKEL PL, PHIL. TRUSSHEAD #6128-148N SERVICE SUPPLY CO.	4.00	EA
CSM-039 (S.M. SCREW)	6 X 1/2" SELF DRILLING #410 S.S. HEX HEAD, CORROSION RESISTANT COATING	36.00	EA
CSM-039 (S.M. SCREW)	6 X 1/2" SELF DRILLING #410 S.S. HEX HEAD, CORROSION RESISTANT COATING	15.00	EA
CSM-039 (S.M. SCREW)	6 X 1/2" SELF DRILLING #410 S.S. HEX HEAD, CORROSION RESISTANT COATING	16.00	EA
CSM-039 (S.M. SCREW)	6 X 1/2" SELF DRILLING #410 S.S. HEX HEAD, CORROSION RESISTANT COATING	20.00	EA
CSM-039 (S.M. SCREW)	6 X 1/2" SELF DRILLING #410 S.S. HEX HEAD, CORROSION RESISTANT COATING	14.00	EA
CWS-024 (DRYWALL ...	8 X 2 PARTIAL BOARD SCREW PARTIAL THREAD 8 X 2 PHILIPS HEAD FLAT HEAD WITH NIBS, BLACK PHOS STEEL	44.00	EA
EBR-145 (BREAKER)	10A, 1-P, ALTECH 1DU10L, UL489 ALTECH CORP. CATALOG #1DU10L, DIN MOUNT (AC) CIRCUIT BREAKER, 10A, 1-POLE, 10kA @ 120/240/277VAC, 50/60Hz, L-SERIES; SCREW TERMINALS, GREEN: OFF/TRIP, RED: ON, LINE/LOAD REVERSIBLE, HIGH MAGNETIC TRIP POINT, UL 489 & CE	1.00	EA
EBR-147 (BREAKER)	2A, 1-P, ALTECH 1DU2L, UL489 ALTECH CORP. CATALOG #1DU2L, DIN MOUNT (AC) CIRCUIT BREAKER, 2A, 1-POLE, 10kA @ 120/240/277VAC, 50/60Hz, L-SERIES; SCREW TERMINALS, GREEN: OFF/TRIP, RED: ON, LINE/LOAD REVERSIBLE, HIGH MAGNETIC TRIP POINT, UL 489 & CE	3.00	EA
EBR-163 (BREAKER)	3A, 1-P, ALTECH 1DU3L, UL489 ALTECH 1DU3L, DIN MOUNT (AC) CIRCUIT BREAKER, 3A, 1-POLE, 10kA @ 120/240/277VAC, 50/60Hz, L-SERIES; SCREW TERMINALS, GREEN: OFF/TRIP, RED: ON, LINE/LOAD REVERSIBLE, HIGH MAGNETIC TRIP POINT, UL 489 & CE	3.00	EA
EBR-169 (BREAKER)	20A, 1-P, ALTECH 1DU20L, UL489 ALTECH 1DU20L, DIN MOUNT (AC) CIRCUIT BREAKER, 20A, 1-POLE, 10kA @ 120/240/277VAC, 50/60Hz, L-SERIES; SCREW TERMINALS, GREEN: OFF/TRIP, RED: ON, LINE/LOAD REVERSIBLE, HIGH MAGNETIC TRIP POINT, UL 489 & CE	1.00	EA
EBU-016 (BUSHING)	SR7N-2 HEYCO INSULATOR BUSHING HEYCO #1244 STRAIN RELIEF BUSHING; BLACK, UL/CSA STRAIGHT THRU-TYPE; FOR USE W/ Ø 0.430" MAX. ROUND CABLES; FOR USE W/ (EPL-030) 7-WIRE CORDSET. REQUIRES A 0.66" WIDE X 0.75" LONG OVAL MOUNTING HOLE.	1.00	EA
EBU-016 (BUSHING)	SR7N-2 HEYCO INSULATOR BUSHING HEYCO #1244 STRAIN RELIEF BUSHING; BLACK, UL/CSA STRAIGHT THRU-TYPE; FOR USE W/ Ø 0.430" MAX. ROUND CABLES; FOR USE W/ (EPL-030) 7-WIRE CORDSET. REQUIRES A 0.66" WIDE X 0.75" LONG OVAL MOUNTING HOLE.	2.00	EA
ECL-001 (CLAMP)	1/4" NYLON CABLE CLAMP	2.00	EA
ECL-031 (CLAMP)	3225T63 304 SS LOOP CLAMP 3/8" DIA. SILICONE RUBBER CUSHIONED, 10/PKG.	2.00	EA
ECL-032 (CLAMP)	5/8" NYLON CABLE CLAMP 5/8" NYLON CABLE CLAMP .05 THICK, .203" HOLE	2.00	EA
ECO-061 (CONNECTOR)	WAGO 231-304/026-000 FOUR CIRCUIT WAGO #231-304/026-000, 4-POLE ORANGE FEMALE CONNECTOR W/ CAGE CLAMP, FOR USE W/ AWG #28-12 WIRE, 300 V/15 A MAX., UL & CSA APPROVED.	3.00	EA
ECO-119 (CONNECTOR)	WAGO 2 CONDUCTORS LEVER CONNECTOR WAGO 221-412, 24-12 GA. WIRE, 600V, 20A, 2 COMMON TERMINALS, UL, CE, EN, 100/BOX	6.00	EA
ECO-120 (CONNECTOR)	WAGO 3 CONDUCTORS LEVER CONNECTOR WAGO 221-413, 24-12 GA. WIRE, 600V, 20A, 3 COMMON TERMINALS, UL, CE, EN, 50/BOX	4.00	EA
ECO-121 (CONNECTOR)	WAGO 5 CONDUCTORS LEVER CONNECTOR WAGO 221-415, 24-12 GA. WIRE, 600V, 20A, 3 COMMON TERMINALS, UL, CE, EN, 25/BOX	1.00	EA
EFS-001 (FLOAT SWI...	GEMS MODEL LS-3, PART #201540, N.O. GEMS MODEL LS-3 FLOAT SWITCH, PART NUMBER 201540, CONFIGURED NORMALLY OPEN (REVERSE FLOAT POSITION FOR N.C. OPERATION); POLYPROPYLENE STEM, MOUNTING AND FLOAT MATERIAL, KYNAR RETAINING CLIP, 3/4" FLOAT DIA., 7/16" ACTUATION, 0.95 MIN LIQ. SPEC. GRAVITY, -40°C TO 100°C, 1/8" NPT MOUNTING, WITH LEAD WIRES	1.00	EA
EHE-030 (HEATER)	175W WATER PAN HEATER, 230 VOL IMMERSION HEATER, REF. DWG. 71-1F-300-51A	1.00	EA
EHE-046 (HEATER)	250W @ 230V U-TUBE HEATER, RED	1.00	EA

Partial Parts List For Job #: 21402-0000

Date Printed: March 21, 2024

For Sales Order: Item: E41L2C8TSWH3

Part #	Description	QTY	UOM
EHE-053 (HEATER)	<p>OPEN AIR HEATER, 250W @ 230V, MARKED WITH RED PAINT, WATLOW PART RBN36B1H-2, 40" OAL X 0.315" DIA., FORMED IN FIGURE 1, A = 17.5625", R = 1", W/ EPOXY SEAL AND 36" SILICONE WIRE LEADS W/ SHRINK BOOT EPOXY SEAL, QUOTE # 710091 See drawing number 2006-1-300-7</p> <p>350W @ 230V U-TUBE HEATER, BLUE</p>	1.00	EA
EHU-073 (HUMIDIFY)	<p>OPEN AIR HEATER, 350W @ 230V, WATLOW PART # RBN36B47H-3, 40" OAL X 0.315" DIA., MARKED W/ BLUE PAINT, FORMED IN FIGURE 1, A = 17.5625", R = 1", W/ EPOXY SEAL AND 36" SILICONE WIRE LEADS W/ SHRINK BOOT EPOXY SEAL, QUOTE # 811185 See drawing number 2006-1-300-7</p> <p>HMP60-A11A1A3A0 INTERCAP HUM. SENS</p>	1.00	EA
EIN-002 (INSULINER)	<p>HUMIDITY SENSOR, VAISALA #HMP60-A11A1A3A0 INTERCAP HUMIDITY SENSOR, MEASUREMENT RANGE: 0-100%RH, 0-1VDC OUTPUT SIGNAL, 3m CABLE, -40 to +60°C OPERATING TEMP, 5...28VDC OPERATING VOLTAGE, ±3% ACCURACY AT 0-40°C</p> <p>422 T&B 1/2", BLUE OR BLACK</p>	2.00	EA
ELD-033	<p>THOMAS & BETTS #422; 1/2" TRADE SIZE, BLUE OR BLACK, INSULINER; UL-APPROVED (UL FILE #E23018); FOR USE W/ (ECO-010).</p> <p>DALI TO PWM CONVERTER (ONE DALI)</p>	1.00	EA
ELD-066	<p>MEAN WELL DAP-04</p> <p>12 X 24 WHITE LED TILE W/COATING</p>	6.00	EA
ELD-071	<p>12" x 24" LED tile with conformal coating. With Samsung LM301B Mid Power LED, 4000K with CR190 SH flux bin. 12 STRINGS OF LEDs WITH 13 LEDs PER STRING (156).</p> <p>LED DRIVER MEANWELL HLG-185H-C1400B</p>	2.00	EA
EPC-032 (USB CABLE)	<p>MEANWELL HLG-185H-C1400B LED DRIVER, CONSTANT CURRENT, 1400mA OUTPUT CURRENT, 200.2W, 71-143V OUTPUT, 90-305VAC INPUT VOLTAGE, WITH DIMMING, -40 TO 90°C OPERATING TEMP RANGE, RoHS, 8.03" x 2.68" W x 1.53" H, IP65 RATING, 3 IN 1 DIMMING</p> <p>PANEL MOUNT USB CABLE, A-A F/M, 3FT</p>	1.00	EA
EPC-044 (COMPUTER)	<p>STARTECH USBPNLAFAM3</p> <p>SIBO, 10" ANDROID TOUCHSCREEN</p>	1.00	EA
EPC-057 (ETHERNET)	<p>SHENZHEN SIBO INDUSTRIAL & DEVELOPMENT CO, LTD. MODEL Q8919 TOUCHSCREEN HMI WITH:</p> <p>ETHERNET CABLE - 5' CAT5 CABLE BLK</p>	1.00	EA
EPL-006 (PLUG)	<p>10X6-12205, 5 FT., BLACK; UL-APPROVED, RoHS-COMPLIANT, WITHOUT PVC BOOT</p> <p>MALE PLUG W/LEAD 18/3 W/54"LD WHITE</p>	2.00	EA
EPL-006 (PLUG)	<p>EMF CORP. M-3-7#PM0030; 3-WIRE CORDSET, W/ 54" LENGTH AWG #18 STRANDED Cu WIRES & 300 V/105°C MAX., 1/2" PRESTRIPPED WHITE PVC INSULATION; UL- & CSA-RECOGNIZED.</p> <p>MALE PLUG W/LEAD 18/3 W/54"LD WHITE</p>	2.00	EA
EPL-041 (PLUG)	<p>EMF CORP. M-3-7#PM0030; 3-WIRE CORDSET, W/ 54" LENGTH AWG #18 STRANDED Cu WIRES & 300 V/105°C MAX., 1/2" PRESTRIPPED WHITE PVC INSULATION; UL- & CSA-RECOGNIZED.</p> <p>9 TERM., AMP MATE-N-LOK #1-480706-0</p>	3.00	EA
EPL-059 (PLUG)	<p>REFERENCE NEWARK #46F2057. TYCO UNIVERSAL MATE-N-LOK #1-480706-0; 9-TERMINAL PLUG, 19 A @ 600 V MAX., 105°C MAX.; USE W/ (ERE-046); UL- &</p> <p>TYCO 172167-1 MINI M-N-L , 4 SOCKET</p>	1.00	EA
EPS-033 (POWER SU...)	<p>172167-1 MINI MATE-N-LOK Connector, Soft Shell; 4; 600 VAC; Plug; Double-Row; Nylon; Plug; Natural</p> <p>100/240V IN, 12VDC OUT, 96W</p>	1.00	EA
ERE-014 (RECEPTACLE)	<p>WAGO 787-1721</p> <p>FEMALE PLUG W/LEAD, 18/3W/12"LD</p>	2.00	EA
ERE-014 (RECEPTACLE)	<p>EMF CORP. #PM0020; 3-WIRE FEMALE RECEPTACLE, W/ THREE (3) 12" LENGTH X AWG #18 STRANDED Cu WIRES & 300 V/105°C MAX., 1/2" PRESTRIPPED WHITE PVC INSULATION; REQUIRES A 1" X 1" SQUARE MOUNTING HOLE; FITS A 0.060" MAX. PANEL THICKNESS; UL- & CSA-RECOGNIZED.</p> <p>FEMALE PLUG W/LEAD, 18/3W/12"LD</p>	2.00	EA
ERL-030 (RELAY)	<p>EMF CORP. PM0020; 3-WIRE FEMALE RECEPTACLE, W/ THREE (3) 12" LENGTH X AWG #18 STRANDED Cu WIRES & 300 V/105°C MAX., 1/2" PRESTRIPPED WHITE PVC INSULATION; REQUIRES A 1" X 1" SQUARE MOUNTING HOLE; FITS A 0.060" MAX. PANEL THICKNESS; UL- & CSA-RECOGNIZED.</p> <p>DPDT,30A.,240VAC</p>	1.00	EA
ERL-060 (RELAY)	<p>S.S., N.O., 3-32 VDC INPUT, 5A, DIN</p>	6.00	EA
ESW-130 (SWITCH)	<p>CARLO GAVAZZI RP1A23D5M1; N.O., 3-32 VDC INPUT, 24-240 VAC OUTPUT, 5 A MAX. LOAD CURRENT, ZERO-SWITCH TRIAC; DIN MOUNT, PUSH-TO-RELEASE WIRING CONNECTIONS; UL- & CSA-APPROVED.</p> <p>DPST ROCKER SWITCH, 20A, 125V</p>	1.00	EA
ETC-052 (TEMP. CON...)	<p>PRCDA1-20F-BB08W</p> <p>DPST ROCKER SWITCH, ON-OFF, 20A @ 125V, 10A @ 250V, PANEL MOUNT, SNAP-IN, BLACK, I-O MARKINGS, RoHS COMPLIANT,</p> <p>3 WIRE RTD, 100 OHMS, 0.00385 PLAT.</p>	2.00	EA
ETC-108 (TEMP. CON...)	<p>WATLOW P/N RB1093901, 128" LEADS</p> <p>INTELLUS CONTROL OVERLAY - CHINA</p>	1.00	EA
ETC-109 (TEMP. CON...)	<p>INTELLUS CONTROL OVERLAY WITH GREEN BACKGROUND FOR CHINA</p> <p>INTELLUSULTRA GEN2, LOW END VER2</p>	1.00	EA
ETE-017 (TERMINAL)	<p>Watlow PSIB-ASIC-0005</p> <p>WAGO 2004-1401, GRAY, 4 TERMINAL</p>	12.00	EA
ETE-018 (TERMINAL)	<p>5 mm TERMINAL BLOCK WIDTH FOR MOUNTING ON DIN RAIL</p> <p>WAGO 2004-1404, BLUE, 4 TERMINAL</p>	5.00	EA
ETE-019 (TERMINAL)	<p>5 mm TERMINAL BLOCK WIDTH FOR MOUNTING ON DIN RAIL</p> <p>WAGO 2004-1407, GROUND, 4 TERMINAL</p>	4.00	EA
ETE-020 (TERMINAL)	<p>5 mm TERMINAL BLOCK WIDTH FOR MOUNTING ON DIN RAIL</p> <p>WAGO 2004-1491, END PLATE</p>	4.00	EA
ETE-021 (TERMINAL)	<p>5 mm TERMINAL BLOCK WIDTH FOR MOUNTING ON DIN RAIL</p> <p>WAGO 2004-402; GRAY TERM. JUMPER</p>	6.00	EA
ETE-022 (TERMINAL)	<p>WAGO #2004-402; LIGHT GRAY, 2-PIN PUSH-IN JUMPER BAR; FOR USE W/ 5 mm WIDE TERMINAL BLOCKS; 32 A MAX.</p> <p>WAGO 249-116, END TERMINAL ANCHOR</p>	12.00	EA
ETE-030 (TERMINAL)	<p>FOR 5 mm TERMINAL BLOCK WIDTH</p> <p>16-CONDUCTOR, 4-TERMINAL SNAP-IN</p>	1.00	EA
EWR-002 (TYRAP)	<p>WAGO #862-0594, BLACK, NO LABELS, MTG REQUIRES TWO Ø5.2mm HOLES SPACED 12mm, PACKS OF 200; UL- & CUL-APPROVED.</p> <p>CABLE TIES, RP14897</p>	2.00	EA
HBS-001 (SPACER)	<p>FEMALE THREADED CIRCUIT BOARD POST</p>	6.00	EA
HHA-008 (HANDLE)	<p>MICRO PLASTICS INC #27FTP00250</p> <p>BLACK SOUTHCO P1-30-101-11</p>	1.00	EA
HHI-029 (HINGE)	<p>ABS LOW GLOSS</p> <p>SOUTHCO EH-5A-4G4-50</p>	2.00	EA
HLO-019 (LOCK)	<p>SERIES 40 BLACK POWDER COATED MOUNTING METHOD: M5 FOUR THREADSTUDS AT .984" BETWEEN CENTER.</p> <p>COMPRESSION LOCK (1/4 TURN)</p>	1.00	EA
	<p>040-1-1-0-42-59, ELECTRAUSA zinc with slotted lock key</p>		

Partial Parts List For Job #: 21402-0000

Date Printed: March 21, 2024

For Sales Order: Item: E41L2C8TSWH3

Part #	Description	QTY	UOM
HMA-010	S-MAGNET 1" X 1/2" X 1/8" THICK BC82CS-S, 1" X 1/2" X 1/8" THICK WITH COUNTERSUNK HOLE FOR #6 SCREW, RECTANGULAR BLOCK MAGNET, N42, NICKEL PLATING, MAGNETIZED WITH THE SOUTH POLE ON THE COUNTERSUNK SIDE	1.00	EA
HNU-045 (NUT)	10-24UNC #18-8 S.S. LOCKNUT SIMILAR TO McMASTER-CARR	3.00	EA
HNU-045 (NUT)	10-24UNC #18-8 S.S. LOCKNUT SIMILAR TO McMASTER-CARR	2.00	EA
HNU-052 (NUT)	M5, .8 COARSE THREAD NUT	8.00	EA
HNU-076 (RIVET NUT)	6-32 RUBBER-INSULATED RIVET NUT McMaster-Carr 93495A120, Rubber-Insulated Brass Rivet Nut, 6-32 Thread, 1/64" - 5/32" Material Thick Range	4.00	EA
HST-016 (PORT PLUG)	1 1/4" PLUG FOR ACCESS PORT	1.00	EA
HWA-001 (WASHER)	NO. 10 PLAIN ZINC	10.00	EA
HWA-012 (WASHER)	1/4" LOCK WASHER ZINC PLATED FASTENAL 1133618	10.00	EA
LST-007 (STAPLES)	15/16 CROWN X 3/4 LEG 16GA. STAPL ROOFING STAPLE	144.00	EA
LST-008 (STAPLES)	15/16 CROWN X 1 LEG 16GA. STAPLE INTERCHANGE WC GALVANIZED STAPLE	278.00	EA
RAC-002 (SUCTION A...	3701 REFRIG. RESEARCH, 5/8 ODF REFRIGERATION RESEARCH #3701, 5/8" ODF SOLDER CONNECTIONS, 2.1 lb R-22 HOLDING CAPACITY	1.00	EA
RBR-001 (RETURN BEND)	1/2" O.D. X .028" W X 1-3/8" #27826 WROT COPPER, MALE ENDS	3.00	EA
RCO-084 (UNIT COOLER)	BOHN TA23BG W/ 230/1/60 FANS BOHN, 230/1/60 Hz, 1.2 AMPS	1.00	EA
RCU-057 (CONDENSIN...	M7GH-H072-CFZ-328, 220/1/50 COMPRESSOR: RRT81C1E-PFZ-150	1.00	EA
RDR-004 (DRIER)	SPORLAN C-083-S DRIER 3/8"	1.00	EA
RFN-001 (FLARE NUT)	FREY #14FS4, 1/4" SHORT TYPE FORGED	1.00	EA
RFN-002 (FLARE NUT)	3/8" SPODEK NS4-6, SHORT TYPE 3/8"	2.00	EA
RFN-003 (FLARE NUT)	NS4-8 1/2" SAE FLARE	3.00	EA
RLI-002 (LIQ. INDICATOR)	SPORLAN #SA-13S; 3/8" ODF SOLDER SPORLAN #SA-13S CATCH-ALL FILTER/DRIER; W/ 3/8 ODF SOLDER CONNECTIONS; 5.2 TON R-22 CAPACITY.	1.00	EA
RRC-028 (RECORDER)	QUENCHARC 0.1MF 100 OHM 72K6759	1.00	EA
RRG-002 (BYPASS VA...	ADRP-3-0/30, 5/8" ODF DISCHARGE BYPASS VALVE	1.00	EA
RTE-004 (TEE)	5/8OD X 1/2OD W-4007 REDUCING TEE MUELLER #W-4007, 5/8" O.D. RUN X 1/2" O.D. BRANCH	1.00	EA
RTE-015 (TEE)	REDUCING, 3/8" R X 1/4" B, #W-40324 MUELLER BRASS #W-40324 REDUCING TEE; 3/8" O.D. RUN X 1/4" O.D. BRANCH.	1.00	EA
RTE-019 (TEE)	TEE, WROT 1/4"NOM SIZE, W-4000 (3/8" OD COPPER TUBE)	1.00	EA
RUN-013 (UNION)	1/4"OD ACCESS FITTING U CAP, 22790 1/4" m. flare access body, solder end has 1/4" ID and 3/8", 7/16" OD steps, 6/PK	1.00	EA
RVH-002 (HAND VALVE)	A-14850 MUELLER 3/8" ODF	1.00	EA
RVQ-001 (BODY)	Q BODY 2 X 4 SAE 168011, 1/4 SAE INLET, 1/2 SAE OUTLET	1.00	EA
RVQ-001 (BODY)	Q BODY 2 X 4 SAE 168011, 1/4 SAE INLET, 1/2 SAE OUTLET	2.00	EA
RVQ-101 (CARTRIDGE)	QC-1 YELLOW Q CARTRIDGE/ID TAG 168072	1.00	EA
RVQ-101 (CARTRIDGE)	QC-1 YELLOW Q CARTRIDGE/ID TAG 168072	1.00	EA
RVQ-102 (CARTRIDGE)	QC-2 GREEN Q CARTRIDGE/ID TAG 168074	1.00	EA
RVQ-200 (ELEMENT)	KT-43-JC 30' ELEMENT 180314	2.00	EA
RVQ-204 (ELEMENT)	KT-43-L3 5' ELEMENT 180331	1.00	EA
RVS-020 (VALVE)	SOLENOID VALVE U141015-220/50 KIP #U141015; 220 VAC/50 Hz, 6 W, 100 MOPD, W/ 1/8" FNPT CONNECTIONS.	1.00	EA
RVS-027 (VALVE)	FEMALE DIN CONNECTOR ASSY GOES W/RVS-028 ON RVS-030 & RVS-031	3.00	EA
RVS-028 (VALVE)	OMKC-1E 220-240/50 COIL GOES W/RVS-027 ON RVS-030 & RVS-031	3.00	EA
RVS-030 (VALVE)	SOLENOID VALVE #XWG, LESS COIL SPORLAN #XWG; W/ EXTENDED 3/8" ODF SOLDER CONNECTIONS	2.00	EA
RVS-031 (VALVE)	SOLENOID VALVE #XUP, LESS COIL SPORLAN #XUP; W/ EXTENDED 3/8" ODF SOLDER CONNECTIONS	1.00	EA
SCL-005 (CLIP)	#66 STAINLESS STL CLIP 0066-000008	16.00	EA
SPI-016 (PILASTER)	S.S. 54" LENGTH, 10060009054 Stainless Steel, 54" long	4.00	EA
SWS-115 (SHELF)	SHELF, 36 1/4" X 27" W/WHITE EPOXY-POLYESTER HYBRID FINISH, DWG. SWS-115(SOLIDWORKS), 41 SERIES	2.00	EA
XBU-001 (BUMPER)	0.78" O.D. X 0.20" THICK, CLEAR	4.00	EA

Partial Parts List For Job #: 21402-0000

Date Printed: March 21, 2024

For Sales Order: Item: E41L2C8TSWH3

Part #	Description	QTY	UOM
XGA-051 (GASKET)	BUMPER SPECIALTIES #BS-41, 0.78" O.D. X 0.20" THICK, CLEAR, SELF-ADHESIVE BACK 39" X 59 3/4" OD, 4 SIDED BLACK	1.00	EA
XGR-001 (GROMMET)	02-126 BLACK MAGNETIC DOOR GASKET, 4 SIDED 1-1/8" OD, 9/32" ID, FOR 7/8" HOLE	2.00	EA
XGR-002 (GROMMET)	3/4" I.D. X 7/8" HOLE	6.00	EA
XGR-002 (GROMMET)	3/4" I.D. X 7/8" HOLE	2.00	EA
XGR-003 (GROMMET)	LAVELLE INDUSTRIES #1022 1-1/4 ID X 1-1/2 HOLE LAVELLE #1047	2.00	EA
XGR-003 (GROMMET)	LAVELLE #1047, 1-1/4" I.D. X FOR 1-1/2" DIA. HOLE X 3/8" THICKNESS X 1/16" GROOVE WIDTH, BLACK 1-1/4 ID X 1-1/2 HOLE LAVELLE #1047	2.00	EA
XGR-003 (GROMMET)	LAVELLE #1047, 1-7/8" O.D. X 1-1/4" I.D. X Ø 1-1/2" PANEL HOLE X 3/8" O.A. THICKNESS X 1/16" GROOVE WIDTH, BLACK RUBBER.	1.00	EA
XGR-003 (GROMMET)	1-1/4 ID X 1-1/2 HOLE LAVELLE #1047	1.00	EA
XGR-003 (GROMMET)	LAVELLE #1047, 1-7/8" O.D. X 1-1/4" I.D. X Ø 1-1/2" PANEL HOLE X 3/8" O.A. THICKNESS X 1/16" GROOVE WIDTH, BLACK RUBBER.	1.00	EA
XGR-005 (GROMMET)	3/8" I.D. X 1/2" HOLE, BLACK VINYL	2.00	EA
XGR-005 (GROMMET)	BLACK VINYL, G1089 PLASTIC	1.00	EA
XGR-005 (GROMMET)	3/8" I.D. X 1/2" HOLE, BLACK VINYL	1.00	EA
XGR-006 (GROMMET)	BLACK VINYL, G1089 PLASTIC	2.00	EA
XGR-006 (GROMMET)	1/2" ID X 13/16" O.D. X 5/8" HOLE	2.00	EA
XGR-013 (GROMMET)	BLACK VINYL, 1/2" I.D. X 13/16" O.D. X 5/8" GROOVE DIA. X 1/16" GROOVE WIDTH X 9/32" THICK PUSH-IN GROMMET 1/4" ID, 5/8" OD	2.00	EA
XGR-022 (GROMMET)	Push-in Grommet 1/4" ID, 5/8" OD WHT LCKING GROMMET FOR 3/8 DIA HOLE	6.00	EA
XLA-013 (LABEL)	McMaster-Carr 5269T232, Locking Grommet for 3/8" Hole Diameter, 3/16" ID, White WAGO 44 X 99 MM, SILVER LABEL	1.00	EA
XLA-016 (DECAL)	WAGO 210-804, 44 MM X 99 MM, 500 MARKERS PER ROLL, WEIGHS 400 GRAMS PER ROLL, SILVER-COLORED PERCIVAL DECAL - LARGE, 19" X 4"	1.00	EA
YDR-002	ACCESS PORT SLEEVE, 1 1/4" ID	1.00	EA
YEL-023 (ELBOW)	COMPONENT HARDWARE R90-1002 Flexible Insulating Conduit Sleeve 90° ELBOW, 1/4 X 1/4 O.D. PUSH-IN	1.00	EA
YEL-026 (ELBOW)	FAST FITTINGS, PART #20-152P, 90° UNION ELBOW, 1/4" X 1/4" PLASTIC PUSH-IN TUBE FITTINGS, 1-140°C OPERATING TEMP RANGE, 0-150psi OPERATING PRESSURE RANGE 90° ELBOW, 1/4 PUSH-IN X 1/8NPTF	1.00	EA
YHI-002 (HINGE)	FAST FITTINGS, PART #20-068P, 90° FIXED ELBOW, 1/4 O.D. PLASTIC PUSH-IN TUBE FITTING X 1/8NPTF, FIXED ELBOW, 1-140°C OPERATING TEMP RANGE, 0-150psi OPERATING PRESSURE RANGE IN LINE KNUCKLE STYLE HINGE	5.00	EA
YHP-012 (HOLE PLUG)	SOUTHCO 96-10-400-11 1/2" WHITE DP-500 PLASTIC 2644	1.00	EA
YHP-015 (HOLE PLUG)	1.01" ID SQ BLK HOLE PLUG	2.00	EA
YHS-001 (SPACER)	Black Polyethylene Square Plug, Push-In, 1.01" Long x 1.01" Wide Inside, 1/2" height, 50D hardness HINGE SPACER 95-3-50-5A	5.00	EA
YMC-006 (MALE CON...)	30% GLASS FILLED NYLON/BLACK MALE CONNECTOR, 1/4 PUSH X 1/8 NPTF	1.00	EA
YNV-005 (VALVE)	FAST FITTINGS, PART #20-054P, MALE CONNECTOR, 1/4 O.D. PLASTIC PUSH-IN TUBE FITTING X 1/8NPTF, STRAIGHT, 1-60°C OPERATING TEMP RANGE, 0-150psi OPERATING PRESSURE RANGE BALL VALVE, 1/4 X 1/4" PUSH-IN	1.00	EA
YWA-001 (WASHER)	FAST FITTINGS, PART #20-900P, PL X PL SHUT-OFF VALVE, 1/4 X 1/4 O.D. PLASTIC PUSH-IN TUBE FITTING, 1-140°C OPERATING TEMP RANGE, 0-150psi OPERATING PRESSURE RANGE .040 PLASTIC SHOULDER WASHER	48.00	EA
	24SW0006 0.040" Nylon 6/6 Shoulder Washer with Natural Finish, Off White		