

TEST CERTIFICATION

This is to certify that the <u>I-36VL</u> chamber serial number serial Scientific, Inc. for Percival Scientific Hong Kong Comeet the specifications listed below:	
The chamber has met the temperature range test of	2° C to 44° C $\pm 0.5^{\circ}$ C (lights off) * 4° C to 44° C $\pm 0.5^{\circ}$ C (lights on) *
The chamber has met the humidity range test of	40% to 85% ± 10% (lights off) *
* - These temperatures and humidity levels are achievable when ambient conditions are 24°C (75°F) or below. Higher ambient temperatures may affect the operating ranges.	
Testing was performed with a calibrated circular chart recorder.	
Date04/29/2024	