

The manual provided is designed to cover the standard configurations for the chamber you ordered. Any special features or optional items ordered are presented below.

Your chamber is equipped with the following:

- Temperature ranges of 2°C to 44°C ± 0.5°C (lights off) and 10°C to 44°C ± 0.5°C (lights on)
- IntellusUltra Controller Tablet / Touch Screen App
- Ultrasonic humidifier (Stulz SCA-1000), dehumidifier with Vaisala HMP60 relative humidity sensor. *
- Light intensity of 405 μmoles / m² / sec as measured 6" from the light fixture at a temperature of 24°C and a light lifetime of 100 hours
- SciWhite® and PetriClear™ LED tiles *
- Open loop dimmable lighting *
- Condensate line routed out rear wall of chamber *
- Air-cooled condensing unit

* - Please refer to the section titled Associated Diagrams, Schematics and Miscellaneous Information for special installation and/or operation instructions.

Your chamber does not include the following options as described in the manual:

- No remote alarm contacts are provided
- No convenience receptacles are provided
- No humidification or dehumidification system is provided.

Note: All supplemental information supersedes any similar information provided in the standard manual.