

## **TEST CERTIFICATION**

This is to certify that the <u>GR-36L</u> chamber serial number <u>34130.02.24</u> provided by <u>Percival Scientific, Inc.</u> for <u>Percival Scientific Hong Kong Co., Ltd.</u> has successfully tested to meet the specifications listed below:
The chamber has met the temperature range test of $2^{\circ}\text{C}$ to $44^{\circ}\text{C} \pm 0.5^{\circ}\text{C}$ (lights off) * $5^{\circ}\text{C}$ to $44^{\circ}\text{C} \pm 0.5^{\circ}\text{C}$ (lights on) *
3 C to 44 C ± 0.3 C (lights on)
N/A The chamber has met the humidity range test of No humidity system ordered
* - These temperatures and humidity levels are achievable when ambient conditions are 24°C (75°F) or below. Higher ambient temperatures may affect the operating ranges.
Testing was performed with a calibrated circular chart recorder.
Date <u>7/5/2024</u>