

TEST CERTIFICATION

This is to certify that the <u>ARC-6L2-PGC</u> chamber serial number <u>34252.01.24</u> provided by <u>Percival Scientific, Inc.</u> for <u>Percival Scientific Hong Kong Co., Ltd.</u> has successfully tested to meet the specifications listed below:
The chamber has met the temperature range test of $2^{\circ}\text{C to } 44^{\circ}\text{C} \pm 0.5^{\circ}\text{C (lights off)} *$
7° C to 44° C $\pm 0.5^{\circ}$ C (lights on) *
N/A The chamber has met the humidity range test of No humidity system ordered
* - These temperatures and humidity levels are achievable when ambient conditions are 24°C (75°F) or below. Higher ambient temperatures may affect the operating ranges.
Testing was performed with a calibrated circular chart recorder.
Date09/24/2024