

The manual provided is designed to cover the standard configurations for the chamber you ordered. Any special features or optional items ordered are presented below.

## Your chamber is equipped with the following:

- Temperature ranges of  $2^{\circ}$ C to  $44^{\circ}$ C  $\pm 0.5^{\circ}$ C (lights off) and  $7^{\circ}$ C to  $44^{\circ}$ C  $\pm 0.5^{\circ}$ C (lights on)
- IntellusUltra Controller Tablet / Touch Screen App
- Light intensity of 680  $\mu moles$  /  $m^2$  / sec as measured 6" from the light fixture at a temperature of 24°C
- SciWhite® LED tiles \*
- Open loop dimmable lighting \*
- Air-cooled condensing unit

\* - Please refer to the section titled Associated Diagrams, Schematics and Miscellaneous Information for special installation and/or operation instructions.

## Your chamber does not include the following options as described in the manual:

- No remote alarm contacts are provided
- No humidification or dehumidification system is provided.

**Note**: All supplemental information supersedes any similar information provided in the standard manual.