

The manual provided is designed to cover the standard configurations for the chamber you ordered. Any special features or optional items ordered are presented below.

Your chamber is equipped with the following:

- Temperature ranges of 2° C to 44° C $\pm 0.5^{\circ}$ C (lights off) and 7° C to 44° C $\pm 0.5^{\circ}$ C (lights on)
- IntellusUltra Controller Tablet / Touch Screen App
- Ultrasonic humidifier (Stulz SCA-1000), dehumidifier with Vaisala HMP60 relative humidity sensor. *
- Light intensity of 800 µmoles / m² / sec as measured 6" from the light fixture at a temperature of 24°C and a light lifetime of 100 hours
- SciWhite® LED tiles with white and red LEDs *
- Closed loop dimmable lighting *
- Air-cooled condensing unit

Your chamber does not include the following options as described in the manual:

• No remote alarm contacts are provided

Note: All supplemental information supersedes any similar information provided in the standard manual.